



HP M1722A/B

CodeMaster XL+
Defibrillator/Monitor
User's Guide

M1722A/B CodeMaster XL+ Defibrillator/Monitor User's Guide

 **HEWLETT®
PACKARD**
HP Part No. M1722-91908
Printed in USA November 1994

Edition 6
E1194

Notice

The information in this document is subject to change without notice.

Hewlett-Packard makes no warranty of any kind on this material, including but not limited to, the implied warranties of merchantability and fitness for a particular purpose. Hewlett-Packard shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

This document contains proprietary information which is protected by copyright. All rights are reserved. No part of this document may be photocopied, reproduced, or translated to another language without the prior written consent of Hewlett-Packard Company.

Before using the instrument, read this guide and become thoroughly familiar with the contents.

Responsibility of the Manufacturer

Hewlett-Packard only considers itself responsible for any effects on safety, reliability and performance of the equipment if:

assembly operations, extensions, re-adjustments, modifications or repairs are done by persons authorized by Hewlett-Packard, and

the electrical installation of the relevant room complies with the IEC or national requirements, and

the instrument is used according to the instructions for use presented in this manual.

NOTE



As with all electronic equipment, radio frequency interference between the defibrillator and any existing RF transmitting or receiving equipment at the installation site, including electrosurgical equipment, should be evaluated carefully and any limitations noted before the equipment is placed in service. Hewlett-Packard assumes no liability for failures resulting from RF

interference between HP medical electronics and any radio frequency generating equipment.

THIS PRODUCT IS NOT INTENDED FOR HOME USE.

This is to certify that this equipment is in accordance with the Radio Interference Requirements of Directive FTZ 1046/84. The German Bundespost was notified that this equipment was put into circulation, the right to check the device for compliance with the requirements was granted.

If Test and Measurement Equipment is operated with unshielded cables and/or used for measurements in open set-ups, the user has to assure that under these operating conditions the Radio Interference Limits are still met at the border of his premises.

CAUTION













Use of accessories other than those recommended by Hewlett-Packard may compromise product performance.

Printing History

March 1992	Edition 1
June 1992	Edition 2
December 1992	Edition 3
July 1993	Edition 4
October 1994	Edition 5
November 1994	Edition 6

© Copyright 1994, Hewlett-Packard Company

Safety Symbols

	Monitor On (Do not confuse with 1 Joule)
	Off (Standby)
	On/Off
	Ground
	Shock hazard
	Caution - See operating instructions
	Meets IEC type BF leakage current requirements and is defibrillator protected.
	Meets IEC type CF leakage current requirements and is defibrillator protected.
	Equipotential (rear of unit, adjacent to a.c. input)
	Protective earth (ground)

Conventions Used in This Manual

WARNING



Warning statements describe conditions or actions that can result in personal injury or loss of life.

CAUTION



Caution statements describe conditions or actions that can result in damage to the equipment or loss of data.

NOTE



Notes contain additional information on usage.

TEXT

Key

LIGHT

represents the labels that appear on the display.

represents keys on the key panel.

represents lighted indicators on the key panel.

Contents

Getting acquainted

Operating Controls and Indications 1-1

Safety considerations 1-8

AC and DC (Battery) Operation 1-10

Battery Life 1-11

Defibrillating

Defibrillating a Patient 2-1

1. Select Energy 2-2

2. Charge 2-2

3. Shock 2-3

After Using the Defibrillator 2-4

Defibrillating with Alternate Paddle Sets 2-4

Monitoring

Using Leads to Monitor 3-1

Preparing the Leads for Monitoring 3-2

Preparing the Patient 3-4

Monitoring Electrodes 3-5

Monitoring 3-6

Heart Rate Alarms 3-7

Printing the Event Summary Record 3-8

Recording 3-9

Automatic Recordings 3-9

Post Shock Data 3-10

Recorder Errors 3-10

External Monitoring 3-11

Performing Synchronized Cardioversion

Performing Cardioversion 4-1

Monitoring During Cardioversion 4-1

Performing Synchronized Cardioversion 4-2

Pacing (Optional)

Using the Pacer 5-2

Defibrillation During Pacing 5-4

SpO₂ Monitoring (Optional)

SaO₂ and SpO₂ 6-1

Application Notes 6-1

Using SpO₂ to Monitor a Patient 6-2

Apply the Sensor to the Patient 6-3

Troubleshooting Sensor Application 6-5

Connect the Sensor to the CodeMaster XL+ 6-6

Start Monitoring 6-7

SpO₂ Readings 6-7

SpO₂ Alarms 6-8

Activating SpO₂ Alarms 6-8

Deactivating SpO₂ Alarms 6-8

Recorder Output 6-9

Troubleshooting

Troubleshooting 7-1

Troubleshooting the Defibrillator 7-2

Troubleshooting the Pacer 7-5

Troubleshooting SpO₂ 7-6

Performing Diagnostics 7-7

Operational Checks 7-8

Every Shift 7-9

Every Day 7-9

Every Week 7-12

Every Three Months 7-12

Every Six Months 7-12

Maintaining the Defibrillator

Changing the Recorder Paper 8-1

Cleaning the Recorder Printhead 8-2

Maintaining the Battery 8-3

Battery Capacity Check 8-4

Replacing the Battery 8-5

Cleaning Exterior Surfaces 8-6

Cleaning and Sterilizing the Internal Paddles 8-7

Steam Sterilizing the Internal Paddles 8-7

Ethylene Oxide Sterilization 8-8

Supplies 8-8

Installation and Setup

Installation A-1

Line Voltage Settings A-1

Installing and Charging the Battery A-1

Loading the Recorder Paper A-3

Connecting Paddles and Patient Cables A-3

Setup A-7

Specifications A-11

Defibrillator A-11

Monitor A-12

Thermal Array Recorder A-12

Size and Weight A-13

Battery A-13

External Pacemaker (Optional) A-14

SpO₂ (Optional) A-14

Calling for Service A-15

List of figures

The CodeMaster XL+ Defibrillator/Monitor	1-1
Defibrillator Operating Controls	1-2
Recorder Operating Controls	1-3
Monitor Operating Controls	1-4
Pacer Operating Controls	1-5
SpO₂ Operating Controls	1-6
Indicator Lights	1-7
Defibrillator Control Panel	2-1
The Monitor Control Panel	3-6
The Pacer Control Panel	5-1
Positioning of the Light Emitters and Photodetector	6-2
Application of the HP M1190A Reusable Sensor	6-4
Connecting the SpO₂ Sensor to the CodeMaster XL+	6-6
CodeMaster XL+ display with SpO₂ reading	6-7
Changing the Recorder Paper	8-1
Installing the Battery	A-3
Connecting External Paddles, Adhesive Pads, or Internal Paddles	A-5
Connecting a Patient Cable	A-7

List of Tables


Defibrillator Operating Controls	1-2
Recorder Operating Controls	1-3
Monitor Operating Controls	1-4
Pacer Operating Controls	1-5
SpO ₂ Operating Controls	1-6
Indicator Lights	1-7
Audible Indicators	1-8
Switchless Internal Paddles Selection	2-8
Switched Internal Paddles Selection	2-9
Cardiac Monitoring Configurations	3-1
3-Wire Electrode Placement	3-2
Lead Formation	3-3
5-Wire Electrode Placement	3-3
Lead Configurations	3-3
Event Summary Record Information	3-8
Automatic Recordings	3-10
External Monitoring Cables	3-11
Failure Messages	7-1
Defibrillator Messages	7-2
Pacer Messages	7-5
SpO ₂ Messages	7-6
Setup Menu 1 Settings	A-9
Setup Menu 2 Settings	A-10

List of Tables

Getting acquainted

This User's Guide provides operational and basic maintenance instructions for safe use and proper care of the Hewlett-Packard M1722A/B CodeMaster XL+ defibrillator/monitor.

Figure 1-1




The CodeMaster XL+ Defibrillator/Monitor

Operating Controls and Indications

The following figures and tables detail the controls and indications on the CodeMaster XL+ defibrillator/monitor.

Figure 1-2



Defibrillator Operating Controls

Table 1-1

Defibrillator Operating Controls






Control	Description
Energy Select/power control	Turns the instrument power on and off and selects energy level.
 button	Charges defibrillator to energy level set on Energy Select control.
 buttons	Administers shock. Labelled  .
 button	Changes operating mode between immediate shock (normal) mode and synchronized with next R-wave shock (Sync) mode.

Figure 1-3




Recorder Operating Controls

Table 1-2

Recorder Operating Controls

Control	Description
Record	Starts and stops the recorder.
Mark	<ul style="list-style-type: none"> When the recorder is on, pressing Mark will annotate the ECG at that point. If the recorder is not on and the unit is set up to do so, pressing Mark will print an ECG strip. <p>See Appendix A for information about setting up the CodeMaster XL+ for recording when you press Mark (Record on Mark).</p>

Figure 1-4




Monitor Operating Controls

Table 1-3

Monitor Operating Controls

Control	Description
Lead Select	Selects an ECG source to monitor.
ECG Size	Changes displayed ECG size.
HR Alarm	Controls HR Alarms.
Review	Prints an Event Summary record. The message "ES" is printed at the top of the ECG strip when you print the Event Summary record. The recorder must be off to print an Event Summary with this key.
QRS Beeper Volume	Controls volume of QRS beeper.
ECG Output	Provides analog 1V/mV output for external monitoring.

Figure 1-5




Pacer Operating Controls

Table 1-4

Pacer Operating Controls

Control	Description
Pacer On	Turns pacemaker on or off.
Rate	Adjusts pacemaker rate (ppm) up or down.
Start/Stop	Starts and stops pacing.
Output	Adjusts pacemaker output current (mA) up or down.
Mode	Changes between fixed and demand pacing modes.

Figure 1-6



SpO₂ Operating Controls


Table 1-5

SpO₂ Operating Controls

Control	Description
SpO ₂ On/Off	Turns pulse oximeter on or off.
SpO ₂ Alarm	Activates, selects and deactivates SpO ₂ alarms.
SpO ₂ input	Connector for SpO ₂ sensor or sensor adaptor cable.

Getting acquainted
Operating Controls and Indications

Figure 1-7



Indicator Lights

Table 1-6

Indicator Lights

Indicator	Description
Charge Done	Indicates that unit is charged and ready to use. Charge lights on key panel and on Apex paddle.
Paddle Contact Indicator (PCI) on Sternum paddle	Indicates how well paddles contact patient.
Sync Light	Indicates that unit is in synchronized shock mode as opposed to defibrillator mode. Flashes off each time an R-wave is detected.
AC Power	Indicates that unit is plugged in to AC power.
Battery Charge	Indicates that the unit is plugged into AC power and that the battery is being charged.

Table 1-7

Audible Indicators

Indicator	Description
Charge Done tone	Sounds when instrument is charged and ready to deliver a shock. Can be disabled in setup.
Auto disarm tone	Sounds during the last ten seconds of the Charge Done tone. Beeps intermittently until disarmed.
QRS beeper	Sounds whenever an R-wave is detected. Volume controlled by front panel adjustment.
CRT alerts	Three beeps each time a message appears on the screen. Can be disabled in setup.
HR alarms	Sounds if the heart rate is above the higher alarm limit or below the lower alarm limit.
SpO ₂ alarm	Sounds if the SpO ₂ level is above the high SpO ₂ alarm limit or below the low SpO ₂ alarm limit.
Shutdown warning	Alternating pitch sounds for 60 seconds when the system is about to turn off. An alert to plug the unit into AC power.

Safety considerations


The CodeMaster XL+ stores high voltage energy and is capable of delivering up to 360 joules of DC energy to a 50 ohm impedance.

- To remove power from the instrument, you must turn the Energy Select control to Off (Standby). Disconnecting the CodeMaster XL+ from an AC outlet will not remove power because the instrument is battery powered.
- To disarm a charged instrument, use one of three methods:
 - Turn the Energy Select control from an energy level setting to the Monitor On or Off (Standby) position.
 - or

- Place the paddles in their holders and depress both Shock buttons.
or
- Leave the unit charged for 60 seconds and it will automatically disarm.

CAUTION



- Do not leave the instrument turned on when it is not in use and it is not plugged into AC power.
- Do not discharge the defibrillator with the paddles shorted together. To do so can cause burning and pitting of the metal paddle contacts.
- Disconnect any other medical electronic equipment from the patient during defibrillation discharge unless labelled as defibrillator protected ()



WARNING



- **Avoid open paddle discharges. Dangerous high voltage exists on the paddles when the defibrillator is discharged. Contact with this high voltage could cause death or serious injury.**
- **Avoid touching any metal surfaces on the instrument during shock.**
- **Avoid connecting the patient to several devices at once, because leakage current limits can be exceeded.**
- **Never touch the bed, the patient, or any equipment connected to the patient during defibrillation.**
- **Keep the CodeMaster XL+ and the immediate area clean and dry at all times to avoid creating potentially dangerous electrical paths.**
- **Never open the instrument case. Dangerous high voltages can be exposed. Only qualified service personnel can service the instrument.**
- **Do not use the defibrillator in a flammable or oxygen-rich atmosphere. This will cause an explosion hazard.**

- **Do not rely entirely upon heart rate alarms. Rate meters on patients can continue to count the pacemaker rates during cardiac arrest or some arrhythmias. Keep pacemaker patients under close observation.**
 - **Avoid moving a charged defibrillator. If the unit is dropped, it may discharge.**
-

AC and DC (Battery) Operation

The CodeMaster XL+ defibrillator can be operated on AC line or DC battery power. The following is a list of AC and battery operating instructions.

- The battery will charge when the instrument is connected to AC power even if the Energy Select switch is in the Off (Standby) position.
- Battery charging is indicated by the **BATT CHRG** light being on.
- A fully depleted battery will recharge to 90% of full capacity in two hours and 100% capacity in 18 hours. To preserve battery integrity, the battery must be fully recharged each time the battery is depleted.
- A new battery or one that has been stored for an extended period requires 24 hours of charging before use.
- When the unit is not in use, connect it to AC power with the monitor turned off. This is to maintain a full battery charge and to prolong battery life.
- To operate on internal battery power only, disconnect the power cord from the AC outlet.
- A fully charged battery will nominally provide fifty 360 joule charge-shock cycles or 2.5 hours of continuous monitoring (15C-40C).

NOTE



Continuous recording will reduce monitoring time available when you are using the unit on battery power.

CAUTION



When the **LOW BATTERY** message is displayed on the monitor, plug the unit into AC power.

From the time the **LOW BATTERY** message is first displayed to when the battery capacity is fully depleted (instrument shutdown), there is typically enough reserve battery capacity to provide either 30 minutes of monitoring, or five 360 joule charge-shock cycles.

If a battery is defective, there is significantly less monitoring or charge capacity available after the **LOW BATTERY** message appears than if the battery is merely depleted.

A unique audible alarm will sound continuously when there are 60 seconds of battery capacity remaining. The instrument will automatically shut off after 60 seconds.

If the battery has been fully depleted, plugging the instrument into AC will immediately restore full operation.

Frequent battery discharges to the low battery level will degrade battery life.

Battery Life

The sealed lead-acid battery used in the CodeMaster XL+ will provide optimum life when the unit is continually connected to AC power when not in use. The battery operates best when it is fully charged after each use. To fully charge a depleted battery requires 18 hours of continuous charge time. Because it is not always practical to allow a full charge cycle between uses, the CodeMaster XL+ can charge a depleted battery to 90% of its capacity in approximately two hours. However, battery capacity and battery life will be reduced if the battery is not allowed to fully charge after each use. For improved battery life, consider ways to reduce the number of instrument uses between full charge cycles.

When the instrument is not plugged into AC power, some current is drawn from the battery to maintain memory and startup logic. Remove the battery if the instrument is to be stored for extended periods (more than one month) without AC power. Note on the instrument that the battery has been removed. After an extended storage period, test the battery according to the battery capacity check as described in *Maintaining the Battery* on page 8-3.

This battery was selected because it provides optimum performance and battery life over a wide range of operating conditions. The life expectancy of this battery is dependent on many variables, including temperature and usage. Periodically check the battery capacity to determine whether to replace it. The battery capacity check is described in *Maintaining the Battery* on page 8-3.

NOTE



When plugged into AC power, the CodeMaster XL+ will function normally with no battery installed, however the time required to charge the defibrillator will increase.

WARNING



If the CodeMaster XL+ is operated without a battery installed, clearly mark the instrument that it does not have a battery and requires AC power to operate. When a CodeMaster XL+ has no battery installed and is plugged into AC power, the front panel **AC POWER** light will be on and the **BATT CHRG** light will be off.

WARNING




Properly dispose of or recycle depleted batteries according to local regulations. Do not disassemble, puncture or incinerate the disposed batteries.

Defibrillating

This chapter contains information about defibrillating a patient and using different paddle sets for defibrillating.

Figure 2-1



Defibrillator Control Panel

Defibrillating a Patient

The following section describes the three steps to defibrillating a patient:

- 1 Select Energy
- 2 Charge
- 3 Shock

1. Select Energy

- 1 Turn the Energy Select control to the desired energy level. The defibrillator is now on.
- 2 Prepare the paddles by following these steps.
 - a. Remove the paddles from their holders by grasping the handles and lifting them straight up.
 - b. Apply Redux® electrolyte paste to the electrode surface of each paddle.

WARNING



Do not allow paste to accumulate on your hands or on the paddle handles to avoid risk of electrical shock.

- 3 Apply the paddles as described below.
 - a. Place the left (Sternum) paddle near the upper sternum just below the patient's right clavicle.
 - b. Place the right (Apex) paddle on the chest just below and to the left of the patient's left nipple, in the anterior-axillary line.
- 4 Rub the paddles lightly against the skin to increase contact between the patient skin and the paddles. Then keep the paddles still to reduce motion artifact on the monitor.

WARNING



Do not spread paste between the paddle electrodes on the chest. The patient can be burned if the paste forms a path between the electrodes.

- 5 Apply 10 to 12 kg of pressure to the paddles.

2. Charge

Press the Charge button on either the right (Apex) paddle or on the instrument front panel.

Wait for the Charge Done indicators. When the unit is armed, the monitor Delivered Energy display shows the available energy in joules.

If the defibrillator does not charge, refer to Chapter 7, **Troubleshooting**.

Resetting the Selected Energy Level

To increase or decrease the selected energy level after pressing the Charge button, perform the following steps.

- 1 Move the Energy Select control to the new energy level.
- 2 Wait for the Charge Done indicators.

3. Shock

To shock the patient, perform the following steps.

- 1 Briefly adjust paddle pressure and placement to optimize patient contact, as registered on the paddle contact indicator (if supplied).
- 2 Verify that no one is in contact with the patient, the monitoring cable or leads, the bed rails, or another potential current pathway.
- 3 Call out "Clear!" to alert personnel to stand away from the patient.

WARNING



Keep hands clear of the paddle electrode edges. Use your thumbs to depress the Shock buttons on the paddle handles.

- 4 Press and briefly hold both Shock buttons (one on each paddle) simultaneously, to deliver energy to the patient.
- If the defibrillator does not shock, refer to Chapter 7, **Troubleshooting**.

NOTE



If you must disarm the charged defibrillator (if countershock is not needed), turn the Energy Select control to Monitor On. Any stored energy will be discharged internally and the available energy on the display will return to 0.

After Using the Defibrillator

If you want to print an Event Summary now, press **REVIEW**. See *Printing the Event Summary Record* on page 3-8. After you use the defibrillator, perform the following steps to prepare the defibrillator for its next use.

- 1 Turn the Energy Select control to Off (Standby).
- 2 Return the instrument to its storage location, and plug the power cord into an AC power outlet. Verify that the **BATT CHRG** and **AC POWER** lights are on.
- 3 Clean all paddles, controls, and cables. Refer to Chapter 8, **Maintaining the Defibrillator** for cleaning instructions.
- 4 Check that sufficient recorder paper and electrolyte paste or defibrillator pads are available for the next use of the defibrillator.

Defibrillating with Alternate Paddle Sets

The CodeMaster XL+ will defibrillate with several different pads/paddles sets.

- Adult/Pediatric Anterior/Anterior External Paddles
- External Adhesive Pads
- Internal Paddles

Performing Pediatric Defibrillation

The CodeMaster XL+ paddle set comes with pediatric paddles. To use the pediatric paddle set, depress the release latch at the front of the standard external paddle set while pulling forward on the adult paddle surface. This action will remove the adult paddle contact surface and expose the smaller pediatric contact surface.

Refer to the *Defibrillating a Patient* section for defibrillation procedures.

WARNING



The clinician must select an appropriate energy level for the pediatric patient. There is no energy limit lockout for the pediatric paddle set.

Defibrillating through External Adhesive Pads

The CodeMaster XL+ has an external pads adaptor that is optional. This adaptor allows defibrillation through external adhesive pads.

External Pads have the following advantages.

- They allow "hands off" defibrillation.
- They provide good quality monitoring.
- You can perform synchronized cardioversion without using an ECG lead set, while monitoring through the pads.
- If the optional pacer is installed, you can switch between the pacing and defibrillation modes of operation quickly.

WARNING



The defibrillator will deliver defibrillator energy levels to an open pads set. The message PADS OFF appears when there is a poor pads-to-patient contact. Check all patient connections if this message appears. If possible, dry off the patient's chest prior to applying the pads. Avoid spreading gel on the patient's body surface. Excess moisture or pads gel between pads can form a path between pads. This reduces defibrillation effectiveness and can possibly burn the patient.

- 1** Attach the pads adaptor cable (HP M1750A/B) to the paddle connector on the front of the defibrillator.
- 2** Slide the paddle connector lock towards the front of the defibrillator to secure the cable.
- 3** Attach the pads to the patient as instructed on the package.


- 4 Connect the pads to the pads adaptor cable. The pads are correctly connected when the locking ring is twisted, locking the ears of the connector to the adaptor cable.

WARNING




Failure to correctly connect the pads to the adaptor cable can result in a failure to deliver energy to the patient.


If the PADS OFF monitor message is displayed, check all patient connections.

- 5 Select pads as the ECG source by pressing **Lead Select** until PADS appears on the display under the heart rate.
- 6 Set Energy Select control to desired energy.
- 7 Press **Charge**.
- 8 Wait for the Charge Done indicators.
- 9 Press both Shock buttons at once to defibrillate. The Shock  buttons for external pads are on the cable connector housing.

Performing Internal Defibrillation

You can perform internal defibrillation using one of the optional internal defibrillation paddle sets.

The Switchless paddle sets attach to the internal paddles adaptor cable (M1740A/B). The Shock buttons  are on the connector housing of the internal paddle adaptor cable.

The Switched paddle set does not require the internal paddles adaptor cable. This paddle set has a single Shock button  on the right-hand paddle handle.

To defibrillate a patient internally, using the switchless paddles, perform the following steps.

Defibrillating
Defibrillating a Patient

- 1 Attach the internal paddles adaptor cable (HP M1740A/B) to the paddle connector on the front of the defibrillator.
- 2 Slide the paddle connector lock towards the front of the defibrillator to secure the cable.
- 3 Select correct paddle set size from Table 2-1, "Switchless Internal Paddles Selection," on page 2-8.
- 4 Attach the internal paddles set to the internal paddles adaptor cable. The pads are correctly connected when the locking ring is twisted, locking the ears of the connector to the adaptor cable.

WARNING



Failure to correctly connect the internal paddles to the adaptor cable can result in a failure to deliver energy to the patient.

- 5 Set the Energy Select control to the desired energy.

WARNING



If the energy switch is set to a level greater than 50 joules, 50J MAXIMUM will be displayed. If **Charge** is pressed, the unit will only charge to 50 joules.

WARNING



For safety and sterility when using the paddles, do not touch the paddle beyond the fingerguard on the handle.


- 6 Apply the internal paddles.
- 7 Press **Charge** .
- 8 Wait for the Charge Done indicators.
- 9 Press both Shock buttons at once to defibrillate. The Shock  buttons are on the cable connector housing.

Table 2-1

Switchless Internal Paddles Selection

HP Part Number	Description
M1740A/B	Internal paddles adaptor
M1741A	Internal paddle set, 7.5 cm diameter
M1742A	Internal paddle set, 6.0 cm diameter
M1743A	Internal paddle set, 4.5 cm diameter
M1744A	Internal paddle set, 2.8 cm diameter

To defibrillate a patient internally, using the switched paddles, perform the following steps.

- 1 Select correct paddle set size from Table 2-2, "Switched Internal Paddles Selection," on page 2-9.
- 2 Attach the switched internal paddles cable to the paddle connector on the front of the defibrillator.
- 3 Slide the paddle connector lock towards the front of the defibrillator to secure the cable.
- 4 Set the Energy Select control to the desired energy.

NOTE



If the energy switch is set to a level greater than 50 joules, 50J **MAXIMUM** will be displayed. If **Charge** is pressed, the unit will only charge to 50 joules.

WARNING



For safety and sterility when using the paddles, do not touch the paddle beyond the fingerguard on the handle.

- 5 Apply the internal paddles.
- 6 Press **Charge**.

- 7 Wait for the Charge Done indicators.
- 8 Press the Shock button on the right-hand internal paddle handle.

Table 2-2

Switched Internal Paddles Selection

HP Part Number	Description
M1784A	Internal paddle set, 7.5 cm diameter
M1785A	Internal paddle set, 6.0 cm diameter
M1786A	Internal paddle set, 4.5 cm diameter
M1787A	Internal paddle set, 2.8 cm diameter

Defibrillating
Defibrillating a Patient

Monitoring

This chapter contains information about monitoring a patient with the CodeMaster XL+ defibrillator/monitor. This chapter also contains details of patient preparation that apply to the synchronized cardioversion and pacing procedures described later in this manual.

The CodeMaster XL+ can be used for either short term or long-term cardiac monitoring. A fully charged battery pack provides a minimum of two and a half hours of continuous monitoring. The power cord can be connected to AC power for unlimited monitoring periods.

Using Leads to Monitor

The CodeMaster XL+'s monitoring functions can be used for cardiac monitoring, elective cardioversion, and pacing (optional). Table 3-1, "Cardiac Monitoring Configurations," on page 3-1 details the different ECG sources that can be used for cardiac monitoring and monitoring applications for which each is suited.

Table 3-1

Cardiac Monitoring Configurations

For this ECG source...	Use this cable type...	In this monitoring application...	Electrode P/N
LEADS, LEAD I,II,III	3-Wire: <ul style="list-style-type: none">● 6 Pin - M1731A● 8 Pin - M1733A● 12 Pin - M1605A/M1500A	<ul style="list-style-type: none">● ECG Normal.● Synchronized Cardioversion.● Pacing.	14445C
LEADS, LEAD I,II,III,aVR,aVL, aVF, and V	5-Wire: <ul style="list-style-type: none">● 6 Pin - M1732A● 8 Pin - M1734A● 12 Pin - M1625A/M1520A	<ul style="list-style-type: none">● ECG Normal.● Synchronized Cardioversion.● Pacing.	40493D

Table 3-1

Cardiac Monitoring Configurations

For this ECG source...	Use this cable type...	In this monitoring application...	Electrode P/N
PADS	Pads adaptor cable: M1750A/B	<ul style="list-style-type: none"> ● ECG Normal ● Synchronized Cardioversion <p>You can use these pads for pacing, however you must select LEADS as the monitor source during pacing.</p>	M1749A
PADDLES	Standard with instrument M1746A	Emergency ECG Monitoring	N/A

Preparing the Leads for Monitoring

The CodeMaster XL+ can be configured to use either a 3-wire or a 5-wire patient cable. Use setup menu 2 as described in Appendix A to select the patient cable type (3-wire, 5-wire).

3-Wire Patient Cable

Table 3-2, “3-Wire Electrode Placement,” on page 3-2 describes typical lead-wire placement using the 3-wire patient cable. Table 3-3, “Lead Formation,” on page 3-3 shows how the individual leads are formed using the individual leadwires.

Table 3-2

3-Wire Electrode Placement

Electrode	Placement
RA/White	Near right midclavicular line, directly below the clavicle
LA/Black	Near the left midclavicular line, directly below the clavicle
LL/Red	Below the left pectoral muscle on the left midclavicular line

Table 3-3

Lead Formation

Lead	+	-	ref
I	LA	RA	LL
II	LL	RA	LA
III	LL	LA	RA

5-Wire Patient Cable

Table 3-4, "5-Wire Electrode Placement," on page 3-3 describes a typical lead-wire placement using the 5-wire patient cable. Table 3-5, "Lead Configurations," on page 3-3 shows how the individual leads are formed using the individual leadwires.

Table 3-4

5-Wire Electrode Placement

Electrode	Placement
RA/White	Near right midclavicular line, directly below clavicle.
LA/Black	Near left midclavicular line, directly below clavicle.
LL/Red	Below the left pectoral muscle on the left midclavicular line.
RL/Green	Below the right pectoral muscle on the right midclavicular line.
V/Brown	As appropriate for the V lead to be monitored (V1 - V6).

Table 3-5

Lead Configurations

Lead	Leadwire Combinations
I	LA - RA
II	LL - RA
III	LL - LA
aVR	RA - .5 (LA + LL)
aVF	LL - .5 (RA + LA)

Table 3-5

Lead Configurations

Lead	Leadwire Combinations
aVL	LA - .5 (RA + LL)
V	V - 1/3 (RA + LA + LL)

Preparing the Patient

Proper application and placement of electrodes is essential for quality ECG monitoring. Good contact between the electrode and the skin reduces the effects of motion artifact and signal interference.

- 1 If necessary, shave hair from the site to ensure good electrode to skin contact.
- 2 Clean the skin with soap and water or with alcohol, then wipe it dry.

NOTE



You can safely monitor a patient during defibrillation. However, monitoring electrodes can become polarized during defibrillation shock, causing the ECG waveform to briefly disappear from the display. You can reduce this effect by using silver-silver chloride electrodes.

- 3 Attach disposable electrodes, perform the following steps.
 - a. Peel the protective backing from the electrode. Be careful to keep adhesive surface free from electrolyte paste.
 - b. Apply the electrodes firmly to the patient's skin, pressing around the entire edge of the electrode.
 - c. Attach snap-on or clip-on leads, assuring good contact between the electrode and the lead end. Tape the lead wire to the skin to prevent the electrode or lead from loosening.
 - d. Plug the patient cable connector into the ECG input connector that is in the lower front of the defibrillator, behind the carrying handle.

NOTE



Be careful to correctly align the cable plug when connecting the patient ECG leads cable to the defibrillator/monitor. Correctly orient the cable plug key with the defibrillator connector slot. If the ECG leads cable falls off or is incorrectly connected, the message **LEADS OFF** appears on the display.

Monitoring Electrodes

Using Pads

Standard pads (HP1749A) allow you to monitor through the pads for defibrillation and synchronized cardioversion. If you wish to use the pads for pacing however, you must attach separate electrodes for monitoring. To use standard pads, perform the following steps.

- 1 Attach pads as instructed on the pads package.
- 2 Attach the pads adaptor cable to the defibrillator.
- 3 Connect the pads to the pads adaptor cable.
- 4 **PADS** is the only selectable ECG source.

Using Paddles

For an emergency evaluation you can monitor a patient's ECG through the paddles electrodes when leads are not attached to the patient.


WARNING



Do NOT use paddles to monitor the ECG during elective cardioversion procedures when the instrument is in synchronous (SYNC) mode. Refer to Chapter 4, Performing Synchronized Cardioversion, for detailed information on performing elective cardioversion.

Monitoring

Figure 3-1



The Monitor Control Panel

To monitor a patient's ECG with the CodeMaster XL+, perform the following steps.

- 1 Prepare the patient for ECG monitoring.
- 2 Turn the Energy Select control to the Monitor On position.
- 3 Press **Lead Select** to select the ECG source. The selected source appears in the upper right corner of the display. For example, PADS appears on the display when it has been selected.
 - o If the message LEADS OFF or PADS OFF appears on the display, inspect the electrodes, patient cable, leadwires, and associated connections. If the selected ECG source is not connected, a dashed line will replace the normal ECG trace.

- 4 Ensure that the ECG size has been automatically adjusted for optimal size. If you wish to reduce the ECG size, press **▽ ECG Size**. The "gain bar" along the left side of the display represents 1 mV of signal amplitude.

NOTE



Autogain allows an initial quick setup when the instrument is turned on. To remove the instrument from Autogain, press **ECG Size**. You then must adjust the ECG size manually.

Adjust the QRS beeper volume to the desired volume.

Heart Rate Alarms

The CodeMaster XL+ provides three configurable pairs of upper and lower heart rate alarm limits. Each pair of heart rate alarm limits can be defined in setup menu 1 as described in Appendix A. While monitoring, you can select and enable any pair of pre-defined limits using **HR Alarm**.

When HR alarms are inactive, the monitor will display a bell symbol with a "v" through it.

To select a pair of HR alarm limits, press **HR Alarm** until the pair of limits you wish to use are displayed. If you do not press the key again, the displayed HR alarm limits become active and the limits are replaced by the bell symbol.

If the HR alarm limits are violated, the HR alarm limits replace the bell symbol, and the violated limit is highlighted. Pressing **HR Alarm** at this point will turn off the HR alarms.

If the HR Alarms are active and you wish to review the limits, press **HR Alarm**. The currently active pair of HR alarm limits are displayed momentarily.

If the HR alarms are active and you wish to select another pair of limits, press **HR Alarm** until the pair of limits you wish to use are displayed. Pressing **HR Alarm** repeatedly cycles through the three pairs of HR alarm limits and the HR alarms inactive choice.

NOTE



HR alarms are automatically turned off when you press **Charge** .

Printing the Event Summary Record

During defibrillator usage, the monitor stores up to 28 ECG strips of critical information called events. Events include all shocks, heart rate alarms, SpO₂ alarms and mark events. Each event record includes date of event, heart rate, ECG source, and size setting as shown in Table 3-6, "Event Summary Record Information," on page 3-8. The time annotated on the ECG strip is within 8 seconds of the recorded event. The message "ES" is printed at the top of the ECG strip when you print the Event Summary record.

Table 3-6

Event Summary Record Information

Event	Event Summary Description
Shock	Shock#, Delivered energy, peak current, and patient impedance
Heart Rate Alarms violation	Heart Rate alarm limits.
SpO ₂ Alarms violation	SpO ₂ alarm limits.

Mark Mark Marker symbol (▼) annotates strip at point **Mark** was pressed.

- To print the Event Summary on the recorder, press **Review** . The recorder must not be printing to print an Event Summary with this key. After printing an event, you must wait 10 seconds before printing another event.
- To stop printing the Event Summary, press **Review** or **Record** .
- To review the Event Summary later, turn the unit on and press **Review** .

NOTE



The Event Summary record is cleared each time the defibrillator is turned on and a new event occurs. This allows you to turn off the defibrillator and return later to review event information such as code statistics.

Turn the defibrillator off between uses to ensure that Event Summary records are patient specific.

Recording

To print a record of the current ECG and of the monitor status, press

Record.

- The upper line of the ECG strip contains a periodic report of monitor parameters (Date, Time, Heart Rate, ECG Source, ECG Size, and Recorder mode).
- The lower line of the ECG strip records asynchronous events such as Shock delivery or Heart Rate Alarm violations.
- Several graphic symbols are used to annotate events such as Shock, HR Alarms, Mark, or Sync.

A 1 mV, 200 ms calibration pulse can be printed on the ECG strip by pressing both arrows on the **ECG Size** key simultaneously.

Appendix A contains a list of configuration settings in setup menu 2 which affect recorder operation. The recorder can be configured for either monitor or diagnostic ECG bandwidth data. Delayed (6 seconds) or non-delayed operation is also configurable.

Automatic Recordings

In setup menu 2, you can enable or disable any of the following automatic recordings:

- Record on Mark
- Record on Charge
- Record on Shock

• Record on Alarms

The automatic recordings for both delayed and non-delayed recorder modes of operation are defined in Table 3-7, "Automatic Recordings," on page 3-10.

Table 3-7

Automatic Recordings

Event	Delayed mode Pre-event time	Delayed mode Post-event time	Non- Delayed mode Post-event time
Mark pressed	6 seconds	3 seconds	3 seconds
Charge	6 seconds	Until Shock or Disarm event	Until Shock or Disarm event.
Shocking the patient	6 seconds	12 seconds	12 seconds
Alarms violation	6 seconds	3 seconds	6 seconds
Disarm	6 seconds	3 seconds	3 seconds
Test discharge	N/A	3 seconds	3 seconds

Post Shock Data

As described in Appendix A, you can enable or disable the recording of post shock statistics in setup menu 2.

- If Post Shock Data is enabled, the defibrillator will record the shock delivery statistics (Actual Delivered Energy, Patient Impedance, Peak Current).
- If Post Shock Data is disabled, the defibrillator will record the energy to which it was charged as the delivered energy. For example, if the unit was charged to 200J, the delivered energy annotation on the ECG strip would be DEL 200J.

Recorder Errors

The message **CHECK RECORDER** appears if an error occurs while recording. If this message appears, check the recorder paper supply. This message may also appear if the recorder door is open.

External Monitoring

The ECG output provides an analog 1V/mV ECG signal for connection to an external monitor. Compatible external monitoring divider cables are listed in Table 3-8, "External Monitoring Cables," on page 3-11.

Table 3-8

External Monitoring Cables

1000:1 voltage Divider Cable Connector Type	HP Part No.
Six Pin	M1782A
Eight Pin	M14482A
Twelve Pin	M1783A

The ECG output is also compatible with the interface to the HP central station ECG input.

NOTE



Do not use the ECG output to synchronize another defibrillator. (The ECG-In to ECG-Out delay is 35 milliseconds.)

CAUTION



The connection of external equipment may increase leakage currents. Always request that local safety personnel verify that multiple connected equipment comply with local regulatory standards before putting such equipment into service.

Performing Synchronized Cardioversion

This chapter contains information about performing synchronized cardioversion on a patient with the CodeMaster XL+ defibrillator/monitor.

Refer to Chapter 3, **Monitoring** for information on patient preparation.

Performing Cardioversion

Treatment for certain arrhythmias require synchronizing a defibrillator shock with the ECG's R-wave. It is essential that this R-wave is detected to avoid inducing ventricular fibrillation.

Monitoring During Cardioversion

Using an External Monitor

WARNING



Whenever possible, we recommend that you perform synchronized cardioversion procedures while directly monitoring the patient through the defibrillator's pads or leads inputs. If you use an external monitor as the ECG source, ask your biomedical technician to verify that the monitor/ CodeMaster XL+ combination will deliver a synchronized shock within 60 ms of the peak of the R-wave. (See the service manual for test procedure.) Use a 1 mV QRS complex with a QRS width of 40 ms. This performance cannot be guaranteed with all commercially available monitors.

When the patient is already connected to bedside monitoring equipment, there is a cable which plugs into the ECG output jack of the bedside monitor and connects to the CodeMaster XL+ for monitoring. To use an external monitor with the CodeMaster XL+, perform the following steps.

- 1 Select Lead I or Lead II on the CodeMaster XL+.
 - 2 Plug the input end of the cable from the monitoring equipment into the ECG input plug on the CodeMaster XL+.
-

Using Pads for Cardioversion

Synchronized cardioversion can be performed with external pads. PADS can be selected as the ECG source for cardioversion because they are a reliable contact for monitoring. To perform synchronized cardioversion with pads, perform the following steps.

- 1 Attach the pads adaptor cable (M1750A/B) to the defibrillator.
- 2 Attach pads to patient as instructed on the package.
- 3 Connect pads to the pads adaptor cable.

If the PADS OFF monitor message is displayed, check all patient connections.

- 4 Select pads as the ECG source by pressing **Lead Select**.
- 5 Set Energy Select control to desired energy.
- 6 Perform the procedure described in **Performing Cardioversion** on page 4-1.

Performing Synchronized Cardioversion

To start cardioversion, perform the following steps.


- 1 Turn the Energy Select control to Monitor On.
- 2 Select the desired ECG lead by pressing **Lead Select**.
- 3 Press **Sync** once to place the instrument in Sync mode. The message SYNC will appear on the display.

NOTE



If the paddles are selected as the ECG source, the message USE LEADS will appear on the display. Although the instrument will allow synchronized shock in paddles ECG mode, leads mode is recommended. Artifact induced by moving the paddles may resemble an R-wave and trigger defibrillator shock.

Cardioversion can be performed with the instrument in Autogain mode. Always inspect the displayed ECG before delivering the countershock, and verify that an R-wave marker (indicating shock point) appears only with each R-wave. If a marker dot does not appear, or if a marker dot is viewed on the T-wave segment of the ECG, follow these instructions.

- Adjust the ECG size by pressing  until the marker dot appears only with each R-wave.
 - Select a different lead or adjust the electrode placement, if necessary, to improve ECG R-wave quality.
- 4 Select the desired energy level with the Energy Select control.

Verifying Defibrillator Operation


Before the cardioversion procedure is started, verify defibrillator operation by performing the following brief test.

- 1 Ensure that the paddles connector is attached and the paddles are in their holders.
- 2 Ensure that the patient cable is connected to the ECG Input jack in the front of the defibrillator and the leads are attached to the patient.
- 3 Turn the Energy Select control to 100 joules.

WARNING



Keep hands clear of the paddle electrode edges. Use your thumbs to depress the Shock buttons on the paddle handles.

- 4 Press  . Wait for the Charge Done indicators, and for the Delivered Energy display to read 100 joules.
- 5 With the paddles pressed firmly into their holders, press and hold both Shock buttons at once. The defibrillator will shock with the next detected R-wave. After the shock, the ECG strip will indicate whether the 100 joule test shock passed or failed.

If the defibrillator does not charge, refer to Chapter 7, **Troubleshooting**.

Applying the Paddles

- 1 Prepare the paddles by performing the following steps.
 - a. Remove paddles from their holders by grasping the handles and lifting straight up.
 - b. Apply a liberal amount of Redux® paste to the electrode surface on each paddle or use defibrillator pads.

WARNING



Do not allow paste to accumulate on the hands, the paddle handles, or the paddle electrodes on the chest to avoid the risk of electrical shock or burns.

- c. Gently rub the electrode surfaces together to evenly distribute the applied paste.
- 2 Apply the paddles to the chest as follows.
 - a. Place the left (Sternum) paddle to the right of the sternum just below the clavicle.
 - b. Place the right (Apex) paddle on the chest just below and to the left of the left nipple, in the anterior-axillary line.
- 3 Rub the paddles slightly against the skin to increase the paddle-to-patient contact.
- 4 Apply 10 to 12 kg (22 - 25 lb) of pressure per paddle.
- 5 Press the Charge button on either the right (Apex) paddle or on the instrument front panel.
- 6 Wait for the Charge Done indicators.

NOTE



If you must disarm the charged defibrillator (if countershock is not needed), turn the Energy Select control to Monitor On. Any stored energy will be discharged internally and the available energy on the display will return to 0.

Resetting the Selected Energy Level

To increase or decrease the selected energy level after the Charge button has been pressed, move the Energy Select control to the new energy level, and wait for the Charge Done indicators.

Delivering the Synchronized Shock

To perform synchronized cardioversion, perform the following steps.

- 1 Verify again that the ECG waveform is stable, and that a marker dot appears only with each R-wave of the cardiac cycle.
- 2 Briefly adjust paddle pressure and placement to optimize contact, as registered on the paddle contact indicator.
- 3 Depress both Shock buttons (one on each paddle) until shock occurs. The defibrillator will shock with the next detected R-wave.
- 4 If additional countershocks are required, readjust the Energy Select control as necessary, and repeat the synchronized cardioversion procedure.

NOTE



Depending on how the unit has been configured, it will either remain in the synchronized shock mode or it will return to defibrillator mode following a synchronized shock. Refer to Appendix A for instructions on configuring the defibrillator for operation after synchronized cardioversion.

If the defibrillator does not shock, refer to Chapter 7, **Troubleshooting**.

After Using the Defibrillator

After using the defibrillator, perform the following steps to ensure that the defibrillator is ready for use.

- 1 Turn the Energy Select control to Off (Standby).
- 2 Return the instrument to storage.
- 3 Plug the power cord into an AC power outlet. Verify that the **BATT CHRG** and **AC POWER** lights are on.


Performing Synchronized Cardioversion
Performing Cardioversion

- 4** Clean all paddles, controls, and cables as necessary. Refer to Chapter 8, **Maintaining the Defibrillator** for information on cleaning the defibrillator.
- 5** Check that there are adequate recorder paper, electrolyte paste, and defibrillator pads for the next use.

Pacing (Optional)

This chapter contains information about pacing with the CodeMaster XL+ defibrillator/monitor.

Figure 5-1



The Pacer Control Panel

The CodeMaster XL+ with the optional pacer can perform external transcutaneous pacing. The pacing option provides demand (synchronous) and fixed (asynchronous) pacing modes. The patient is connected to the pacer by external adhesive pads. The patient can be paced and defibrillated through the same set of pads.

WARNING



While pacing, avoid touching the gelled area of the pacing pads or the patient to prevent electrical shock.

Use only HP recommended pads (M1749A) with the external pacer option. The CodeMaster XL+ delivers pacer pulses through a low-impedance multifunction pad. The CodeMaster XL+ does not pace effectively with high-impedance, pace-only electrodes.

Do not use pads for more than eight hours of continuous pacing.

The CodeMaster XL+ will pace on battery power alone. However, whenever possible, plug the CodeMaster XL+ into AC power while pacing.

Using the Pacer

To use the pacer, perform the following steps.

- 1 Apply the pads (M1749A) as instructed on the package.
- 2 Attach monitoring electrodes as instructed in **Using Leads to Monitor** on page 3-1.
- 3 Attach the patient cable to the CodeMaster XL+'s ECG Input connector.
- 4 Attach the patient cable leads to the monitoring electrodes.
- 5 Attach the pads adaptor cable (M1750A/B) to the defibrillator output connector. Pull the latch connector toward the front of the defibrillator to lock the connector in place.
- 6 Attach the pads to the pads adaptor cable and turn the twist lock.
- 7 Turn the Energy Select control to the Monitor On position.
 - o If the message **NO PADDLES** is displayed, check that the pads adaptor cable connector is properly seated and latched.
 - o If the message **PADS OFF** is displayed, check the pads connection to the patient and to the pads adaptor cable.
- 8 Press **(Pacer On)** to turn the pacer on. Pacer parameters will now be displayed at the bottom of the display (**PACER STOP, DEMAND MODE, 70 PPM 30 MA**). The rate (ppm) and output (mA) settings for when the pacer is turned on can be selected in setup menu 1. The original rate and output settings from the factory are 70 ppm and 30 mA. The pacer is always in Demand mode when it is turned on.

NOTE



At this point, no pacer pulses are being delivered to the patient. The pacer must be started before the pacer pulses are delivered at the selected rate (ppm) and output (mA).

- 9 Press **Lead Select** to select the best Lead for monitoring while pacing. You can only select Leads as the ECG source when the pacer is on.
 - o If the message **LEADS OFF** is displayed, check all patient cable connections.
- 10 Press **Rate** to adjust the rate. The selected rate (PPM) is displayed on the monitor.
- 11 Press **Mode** to select the pacing mode (Demand Mode/Fixed Mode). The selected mode is displayed on the monitor.
 - o When in Demand mode, the pacer will only deliver pacer pulses when the patients heart rate is lower than the selected pacer rate.
 - o When in Fixed mode, the pacer will deliver pacer pulses at the selected pacer rate.

WARNING



Use Demand pacing mode whenever possible. Use Fixed pacing mode (Asynchronous) for cases when reliably monitoring the patient is impractical. For example, use Fixed mode when there is motion artifact or other ECG noise that makes R-wave detection unreliable.

- 12 Press **Start/Stop** to start pacing. The monitor will now display the message **PACING** as well as the selected mode, rate and output.

The pacer will not start pacing if there is a problem with either the pacing pad connections or the monitoring electrode connections.

If there is a problem with the pacing pad connection, the message **ATTACH PADS** will be displayed briefly when you press **Start/Stop**.

- 13 Verify that the pacer pulses are well positioned in the diastole.

14 Increase output (mA) by pressing **Output ▲** until the beat is captured. Selecting an alternate lead can help you to determine capture.

15 To set the lowest possible output level to capture, decrease the current by decrements of 5 mA by pressing **Output ▼**.

NOTE



If the monitoring ECG lead falls off while pacing in Demand mode, the pacer will stop delivering pulses and the messages **PACER STOP** and **LEADS OFF** will appear. To resume pacing, reattach the lead and press **Start/Stop**. Pacing in Fixed mode does not require leads to be attached for the pacer to deliver pulses.

If a pacing pad comes off during pacing, the pacer will stop delivering pulses and the messages **PACER STOP** and **PADS OFF** will appear. To resume pacing, reattach the pad and press **Start/Stop**.

WARNING



HR meters and HR alarms function during pacing, but they can be unreliable. The HR meter attempts to count QRS activity in both Demand and Fixed pacing modes. Observe the patient closely while pacing. Do not rely on HR alarms or the indicated Heart Rate as a measure of the patient's health.

Defibrillation During Pacing

If the patient must be defibrillated during pacing, perform the following steps.

- 1 Set the desired energy level with Energy Select control.
- 2 Press **Charge**. The defibrillator will automatically turn off the pacer and start charging. The pacer status messages will be cleared and replaced with the defibrillator Delivered Energy Display.
- 3 Wait for the Charge Done Indicators.
- 4 Call out "Clear!" to alert personnel to stand away from the patient.

Pacing (Optional)
Using the Pacer

- 5 Press and briefly hold both Shock buttons. The shock will be delivered through the multifunction pads.**

After the shock, the pacer remains off. Resume pacing if it is required.

Pacing (Optional)
Using the Pacer



SpO₂ Monitoring (Optional)

SpO₂ Monitoring (Optional)

SpO₂ monitoring gives information on both cardiac and respiratory systems, and details of oxygen transportation in the body. It is widely used because it is non-invasive, continuous, easily applied and painless.

The quality of SpO₂ measurements depends on careful application of the sensor. Read the following section, "Application Notes" (page 6-1), to understand the importance of sensor application. For more detailed information, refer to *SpO₂ Concepts* (HP part number M1722-93950).

You can use the SpO₂ monitor with sensors made by other manufacturers as well as with HP sensors. For a list of approved sensors, see the *Sensor Guide* (HP part number M1722-93970).

SaO₂ and SpO₂

Hewlett-Packard is adopting the convention of referring to the SpO₂ parameter. Previously it was referred to as SaO₂.

SaO₂ is the term used to indicate the oxygen saturation of arterial blood.


SpO₂ is the term used to indicate the oxygen saturation of arterial blood *as measured by pulse oximetry*.

Application Notes

The pulse oximetry method used for measuring SpO₂ uses LEDs (light emitting diodes) to transmit red and infrared light through suitable peripheral areas of the body, typically the foot in neonates or the finger in adults. The oxygen saturation is gauged by measuring the "redness" of the blood in the arterial pulse.

A photodetector positioned opposite the light emitter compares light absorption before and after pulsation to provide measurements that are displayed on the monitor. If there is no pulse, measurements cannot be made. See the following figure.

Figure 6-1



Positioning of the Light Emitters and Photodetector

For accurate measurements, the following conditions must apply:

- All transmitted light must pass through the extremity to the detector.
- The patient must have at least a minimum pulse.
- The light emitter and the photodetector must be opposite each other.

Using SpO₂ to Monitor a Patient

There are three types of sensors:

Disposable

Disposable sensors should be used once only and then discarded. However, they can be relocated to a different application site on the patient if the first location does not give the desired results. Disposable sensors must not be reused on different patients.

Semi-Disposable

Semi-disposable sensors can be reused, but the adhesive wrap must be discarded after each use. Semi-disposable sensors are recommended for single-patient use only.

Reusable

Reusable sensors can be reused on different patients.

Before you start SpO₂ monitoring:


- 1 Estimate the patient's weight, and determine the best site for the sensor.
- 2 Use the *Sensor Guide* to select the correct type and size of sensor for the identified location.
- 3 Prepare the sensor:
 - Disposable: Remove protective backing.
 - Semi-Disposable: Apply a new adhesive wrap to the sensor.
- 4 Apply the sensor to the identified location.
- 5 Connect the sensor to the monitor. To connect sensors from other manufacturers you need the HP M1900B Connector Cable.

Apply the Sensor to the Patient

HP supports the use of many sensors. Use the *Sensor Guide* to find the sensor which is best for your case. Follow the manufacturer's guidelines for applying and using the sensor.

Application of the Reusable Sensor

Figure 6-2



Application of the HP M1190A Reusable Sensor

Push the sensor over the fingertip so the cable lies on the back of the hand, and secure the cable to the wrist with the wrist-strap supplied. Make sure the finger is not pinched in the end of the sensor. This ensures that the light sources in the sensor lie over the base of the fingernail, giving the best measurement results. If the sensor is not in the correct position, inaccurate readings result. In extreme cases, the instrument displays dashes instead of an SpO₂ reading. **When correctly positioned, the end of the finger just touches the end of the sensor.**

CAUTION



When non-HP SpO₂ sensors are used, application must be consistent with the manufacturer's own guidelines.

WARNING



Prolonged, continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering or pressure necrosis, particularly on neonates and on patients with impaired perfusion and varying or immature skin morphology. Specific attention must be given to sensor site inspection for changes in skin quality, proper

optical path alignment and attachment. Check the application site at regular intervals — at least every two hours — and change the site if any compromise in skin quality should occur. More frequent checking may be required due to an individual patient's condition.

Troubleshooting Sensor Application

Failure to apply the sensor properly may cause incorrect measurement of arterial oxygen saturation.

Do not use a damaged sensor or one with exposed electrical circuits.

Patient Movement

Make sure that the application site chosen does not move excessively, which may adversely affect the performance of the sensor. You may have to replace the sensor to ensure good adhesion, or you may have to choose another application site.

Inspecting the Application Site

Inspect the SpO₂ sensor site at least once every 2 hours to ensure adhesion, skin integrity, and correct alignment of the light emitter and photodetector. Should alterations of skin integrity occur, remove the sensor and reapply at another recommended site. Avoid application of the sensor to edematous or fragile tissue. Check circulation distal to the sensor site routinely.

Circulation at Application Site

Wrapping the tape too tightly, or using supplemental tape, can cause venous pulsations that could potentially lead to inaccurate saturation measurements. Therefore, do not wrap the adhesive too tightly and do not use additional tape to secure the sensor. High positive intrathoracic airway pressures, valsalva maneuvers, or other consequences of impaired venous return may also cause venous pulsations.


Only use adhesive wraps recommended by Hewlett-Packard.

Avoid placing the SpO₂ sensor on any extremity with an arterial catheter, blood pressure cuff, or intravascular venous infusion line.

Connect the Sensor to the CodeMaster XL+

When you have applied the sensor to the patient, plug the disposable and semi-disposable sensors into the connector cable and plug this cable into the SpO₂ socket on the lower right of the CodeMaster XL+. Plug the HP M1190A sensor directly into the SpO₂ socket of the CodeMaster XL+. The plug is keyed and is color-coded blue to distinguish it from the white ECG socket.

Figure 6-3




Connecting the SpO₂ Sensor to the CodeMaster XL+

CAUTION




Do not force the SpO₂ connector into the ECG input socket. Doing so may damage the pins on the cable connector.

Start Monitoring

Turn the defibrillator on, if necessary, by turning the Energy Select control to Monitor On. Press the  button to display the SpO₂ reading in the upper right corner of the display.

SpO₂ Readings

Figure 6-4



CodeMaster XL+ display with SpO₂ reading

The **pulse amplitude indicator** shows the quality of the SpO₂ signal. Since it is derived from the patient's plethysmograph signal, it varies with the pulse of the patient. If the patient has a very low signal the pulse amplitude indicator does not vary through its full range. If the signal is noisy, the pulse amplitude indicator does not vary rhythmically with the pulse.

The pulse rate is derived from the pulse oximeter. It should correlate closely with the patient's heart rate.


SpO₂ Alarms

Activating SpO₂ Alarms

There are three preset high/low SpO₂ alarm limits: 100/90, 100/85 and 100/80.

Press the  button repeatedly to cycle through the alarm options

and the no-alarm option. Stop when you see the alarm you would like to


choose and after three seconds that alarm will take effect. A  symbol replaces the limits to show that the alarm is active. To review the limit set,

press the  button.

If the SpO₂ level falls below the low alarm limit, an alarm sounds and the violated limit is highlighted.

NOTE



SpO₂ alarms are automatically turned off when you press  .

Deactivating SpO₂ Alarms

Press the  button. The  symbol to the right of the SpO₂

display indicates that the alarms are off.

Recorder Output

After an alarm event, the recorder prints a strip. The bottom of the strip shows the alarm violation, and the top of the strip shows the SpO₂ reading.

You may print an event summary record, which contains SpO₂ information, as described in “Printing the Event Summary Record” (page 3-8).

Troubleshooting

This chapter contains information about troubleshooting and performing diagnostics on the CodeMaster XL+defibrillator/monitor.

Troubleshooting

This section provides information about messages that appear on the display.

Table 7-1

Failure Messages

Message	Possible Solutions
DEFIB FAILURE	Use a back-up defibrillator and call service.
MONITOR FAILURE	Use a back-up defibrillator and call service. If there is no alternative, try to use the defibrillator without the monitor. Use Charge Done indicators to verify defibrillator functionality.
SYSTEM FAILURE	Use a back-up defibrillator and call service. If there is no alternative, try to use the defibrillator. Use Charge Done indicators to verify defibrillator functionality.

Troubleshooting the Defibrillator

Table 7-2

Defibrillator Messages

Message	Cause	Possible Solutions
NO PADDLES	Paddle set is not connected to the defibrillator.	Attach external paddles, internal paddles, or pads adaptor as required.
LEADS OFF	Leads are not securely attached to the patient or the cable is not connected to the defibrillator.	Attach lead.
USE LEADS	<p>Paddles are connected to the defibrillator and:</p> <ul style="list-style-type: none"> ● Sync was pressed when PADDLES is selected. ● Defibrillator is in sync mode and PADDLES is selected. ● Defibrillator is in defib mode and PADDLES is selected when leads are attached to the patient. ● Defibrillator is in defib mode, PADDLES is the ECG source when attaching leads to the patient. 	<ul style="list-style-type: none"> ● Press Sync to remove the defibrillator from sync mode. ● Select PADS or LEAD I, II, or III instead of PADDLES. ● Select PADS or LEAD I, II, or III instead of PADDLES. ● Select PADS or LEAD I, II, or III instead of PADDLES.
50J MAXIMUM	Internal paddles are attached and the Energy Select control has been turned past 50 joules.	CodeMaster XL+ will charge to 50 joules.
LOW BATTERY	CodeMaster XL+ is not connected to AC power and the battery voltage is below the low battery threshold.	Connect the unit to AC power.
SETUP LOST		Call for service

Table 7-2

Defibrillator Messages

Message	Cause	Possible Solutions
DEFIB DIS-ARMED	CodeMaster XL+ has internally discharged energy and is now disarmed.	
CHECK RECORDER	The recorder door is open. or The recorder is out of paper.	Close the recorder door or replace the paper roll as described in "Changing the Recorder Paper" (page 7-1).

If the Defibrillator does not Charge

If the defibrillator does not charge, perform the following steps.

- 1 Verify the proper setting of the Energy Select control.
- 2 If it is correct, follow these steps.
 - a. Turn the Energy Select control to Off (Standby), and then back to the desired energy setting.
 - b. Press **Charge** again.
- 3 If the instrument remains unable to charge, turn the Energy Select control to Off (Standby) and use a backup defibrillator.
- 4 Press **Review** to print an event summary and keep any ECG strips from the defibrillator for later evaluation.
- 5 Alert appropriate service personnel.

If the Defibrillator does not Deliver a Shock

If the defibrillator does not deliver a shock, perform the following steps.

- 1 Make sure the unit is not in synchronized shock mode. The **SYNC** light is on when the unit is in synchronized shock mode.
- 2 If the unit discharges internally (that is, the energy display decrements slowly then beeps three times and displays the screen message DEFIB

DISARMED) verify proper connections from the patient to the defibrillator. This includes connections to the pads/paddles adaptor cable if one is being used. Also check for worn or broken areas along the cables.

- 3** Press **Charge** again, wait for the Charge Done indicators.
- 4** Press the Shock buttons again.
- 5** If the unit remains unable to shock, turn the Energy Select control to the Monitor On or Off (Standby) position and use a back-up defibrillator.
- 6** Press **Review** to print an event summary and keep any ECG strips from the defibrillator for later evaluation.
- 7** Alert appropriate service personnel as soon as possible.

Troubleshooting the Pacer

Table 7-3


Pacer Messages

Message	Cause	Possible Solutions
PACER FAILURE	The unit detected a delivered current error.	<ol style="list-style-type: none"> 1 Change pads and make sure pads adaptor cable is properly connected. 2 Turn the Energy Select switch to Off (Standby), then turn back to Monitor On. 3 Restart the pacer. If PACER FAILURE happens again, use a back-up pacer and call for service.
PACER OUT-PUT LOW Pacer current high-lighted on the display.	The pacer cannot deliver the required current.	<ol style="list-style-type: none"> 1 Stop the pacer. 2 Change the pads and check the pads adaptor cable for proper connection. 3 Restart the pacer.
STOP PACER	An attempt was made to change pacing mode while pacing.	Stop the pacer before you change the pacing mode.
NO PADS	Pacer is on and one of the following occurs: <ul style="list-style-type: none"> ● No pads adaptor attached. ● Paddles attached. ● Internal paddles attached. 	Attach pads for pacing.
PADS OFF	Pads adaptor cable is attached and a pad is off.	Attach pad for pacing. If pads are already applied to patient, change to a new set of pads.
ATTACH PADS	Pads adaptor cable is attached. (Charge) is pressed and pads are not attached. (Pacer) is pressed and pads are not attached.	Attach the pacer pads.

Troubleshooting SpO₂

Table 7-4

SpO₂ Messages

Message	Cause	Possible Solutions
SPO2 FAILURE	The unit detected a failure in SpO ₂ subsystem hardware. This failure does not affect other parts of the instrument.	Press  twice to power cycle the SpO2 option. If the error happens again, call for service.
SPO2 SENSOR FAIL	The sensor or adaptor cable is broken.	Replace the sensor or adaptor cable.
SPO2 CABLE OFF	Sensor or cable is disconnected.	Check SpO ₂ connections between CodeMaster XL+ and the sensor cable or adaptor cable.
Dashes appear on display instead of SpO ₂ reading.	Can't derive measurement because: <ul style="list-style-type: none"> • The sensor is not on the patient. • No pulse is detected. • The sensor is incorrectly positioned. • The sensor is defective. 	Check the patient for a pulse. Reapply the sensor, and make sure it is correctly positioned. If it doesn't work, replace the sensor.
SPO2 NOISY SIGNAL	Irregular pulse patterns. Patient motion.	Reapply sensor. Consider using a different sensor site.
SPO2 LIGHT INTERF	Too much interference from external light. Damage to sensor or adaptor cable.	Reapply sensor. Turn off lights in the room. If other options do not work, replace the sensor.
SPO2 LOW SIGNAL	Bad connection to patient, or patient has poor perfusion.	Check the patient for poor perfusion. Reapply disposable and semi-disposable sensors. Readjust reusable sensors. Consider using a different sensor site.

Performing Diagnostics

The following procedure allows complete functional inspection of the CodeMaster XL+ defibrillator/monitor.

- 1 Plug the power cord into an AC outlet. Verify that the **(BATT CHRG)** and **(AC POWER)** lights are on.
- 2 Turn the Energy Select control to the Monitor On position. The monitor trace will appear within ten seconds.
- 3 Press **(Lead Select)** until Lead I is displayed and verify that the message **LEADS OFF** appears, indicating that one or more leadwires are not connected.
- 4 Press **(Lead Select)** to return to the paddles ECG selection.
- 5 With the paper installed, press **(Record)** once to turn on the recorder.
 - a. Allow the recorder to run for approximately 20 seconds and check that Date, Time, HR (heart rate), PADDLES (ECG source), and AUTO-GAIN (ECG gain mode) are noted on the ECG strip.
 - b. Press **(Mark)** and verify that the mark (**▼**) symbol is printed. It will be delayed 6 seconds if the recorder is in the delayed mode.
 - c. Press **(Record)** to stop the recorder.
- 6 Press **(HR Alarm)**. Verify that the configured alarm limits appear briefly, then are replaced by the bell symbol bell in the upper right of the display. With no ECG signal, an audible alarm tone should sound within four seconds.
- 7 Press **(HR Alarm)** to turn off the alarms.
- 8 Verify that the adult paddle electrodes are installed.
- 9 Turn the Energy Select control to 100 joules.
- 10 Leaving the paddles in their holders, press either Charge button.

- o The Charge Done indicators should occur within two seconds when operated with a fully charged battery.
- o The Delivered Energy display should indicate 100 joules.

WARNING



Keep hands clear of the paddle electrode edges. Use your thumbs to depress the Shock buttons on the paddle handles.

- 11 Grasp the paddle handles, and without removing the paddles from their holders, press both Shock buttons simultaneously. A brief automatic recorder run prints the delivered energy test report.
- 12 Press **(Sync)** to place the instrument in Sync mode.
- 13 Verify that the messages SYNC and USE LEADS appear on the display.
- 14 Press **(Lead Select)** once to select Lead I, and verify that the message USE LEADS no longer appears. The message LEADS OFF should now appear, indicating that one or more leads are not connected.
- 15 Turn the Energy Select control to Off (Standby).

The defibrillator is ready for use if it passes the above checklist.

A more extensive test of defibrillator/monitor functionality can be performed using the diagnostic service mode described in the *M1722/M1723 CodeMaster XL Series Defibrillator/Monitor Service Manual*.

Operational Checks

These checks are intended to briefly verify the proper operation of the CodeMaster XL+ defibrillator/monitor. Regularly perform a test routine incorporating the following checks along with visual inspection of all cables, paddles, and controls.

Every Shift

Perform the following checks every shift.

- Verify that the instrument is connected to AC power, and that the (BATT CHRG) and (AC POWER) lights are on.
- Check for adequate thermal paper in the recorder.
- Check for ECG leads, electrodes, and adequate REDUX® electrolyte paste or defibrillator pads.

Every Day

- Visually check AC power cord for wear.
- Visually check the patient cables, paddles, cables, and pads adaptor cables for wear, insulation nicks, and other damage.
- Ensure that the (BATT CHRG) and (AC POWER) lights are on. If the unit is plugged in and the (AC POWER) light is not on, the power cord may have a broken wire.
- Check that the (BATT CHRG) and (AC POWER) lights go off when the unit is unplugged.
- Perform the Delivered Energy and Shock Button Functional Test, which follows.
- If your instrument has the external pacing option, perform the Quick Pacer Functionality Test as described in **Quick Pacer Functionality Test** on page 7-11.

Delivered Energy and Shock Button Functional Test

To check the instrument with the paddles, perform the following steps.

- 1 Turn the Energy Select control to the 100 joules position.
- 1 Verify that the adult paddle electrodes are installed.
- 2 Push the paddles completely into their holders (Apex paddle in right pocket, Sternum in left) and press either Charge button. Wait for the Charge Done indicators.

WARNING



Keep hands clear of the paddle electrode edges. Use your thumbs to depress the Shock buttons on the paddle handles.

- 3 With the paddles in their holders, grasp the paddle handles and press the Apex paddle Shock button. Verify that the defibrillator does not discharge.
- 4 Release the Apex paddle Shock button, then press the Sternum paddle Shock button. Verify that the defibrillator does not discharge.
- 5 Press **Sync** to place the defibrillator in sync mode.
- 6 Press and hold both Shock buttons. Verify that the defibrillator does not discharge.
- 7 Press **Sync** again to remove the defibrillator from sync mode.
- 8 With the paddles in their holders, press and briefly hold both Shock buttons at once.
- 9 The recorder will print a test report.

To check the instrument with the pads adaptor cable, perform the following steps.

- 1 Connect the Test Load (M1781A) to the Pads Adaptor Cable (M1750A/B).
- 2 Turn Energy Select to 100 joules.
- 3 Press Charge. Wait for the Charge Done indicators.

- 4 Press each Shock button independently. Verify that the defibrillator does not discharge.
- 5 While depressing Charge, press and briefly hold both Shock buttons.
- 6 The recorder will print a test report.

NOTE



Notify Service Personnel if the ECG strip prints **TEST 100J FAILED** or if any of the Shock button tests fail.

- 7 Disconnect the Test Load (M1781A) from the Pads Adaptor Cable (M1750A/B).

Quick Pacer Functionality Test

- 1 Connect the pads adaptor cable (M1750A/B) to the defibrillator and the test load (M1781A) to the adaptor cable. Turn the unit on by turning the Energy Select switch to Monitor On.
- 2 Press **Pacer On** (ignore the **LEADS OFF** message if it occurs).
- 3 Put the pacer into fixed mode by pressing **Mode**.
- 4 Adjust the pacer current to 30mA by pressing **Output**. Adjust the pacer rate to 60ppm by pressing **Rate**.
- 5 Start the pacer by pressing **Start/Stop**, and start the recorder by pressing **Record**.
- 6 Verify that the pacer pulses are shown on the recorder strip approximately every five large boxes (if the recorder is in delay mode it will take several seconds before the pacer pulses appear). Allow the pacer to run for 10-12 seconds.
- 7 Turn the pacer off by pressing **Pacer On** and stop the recorder by pressing **Record**. Disconnect the test load (M1781A) from the adaptor cable.

Notify service personnel if:

- **PACER FAILURE** is displayed on the monitor.

- The unit beeps three times and displays **PACER OUTPUT LOW**.
- The pacer pulses are not shown on the recorder strip as described in the test.

Every Week

Perform the following checks on the internal paddle set every week.

- Check for excessive residue from sterilization on the paddle set and clean as needed. Oxidation can be an indication the paddle set is old and must be replaced.
- Check for pitting or discoloration on the electrode surfaces. Polish or replace as required.
- Ensure that the cable, connector, and electrodes have no cracks in the insulation.

Every Three Months

Have the cable set tested for electrical continuity every three months.

Every Six Months


See the *M1722/M1723 CodeMaster XL Series Service Manual* for extensive electrical, operational and safety tests to be performed by a qualified Biomedical Equipment Technician (BMET) or equivalent service technician every three to six months.

Maintaining the Defibrillator

This chapter contains information about maintaining and cleaning the CodeMaster XL+ defibrillator/monitor.

Changing the Recorder Paper

Figure 8-1



Changing the Recorder Paper

To change the recorder paper, perform the following steps.

- 1 Slide the recorder door to the right of the defibrillator as shown in Figure 8-1, on page 8-1. The paper platen will tilt up.
- 2 Pull up on the plastic removal tag to remove the empty or low paper roll.

- 3 Place a new roll of thermal paper (HP 40457C/D) in the recorder so that the paper unrolls from the top of the roll and the grid faces down as it comes out of the recorder.**
- 4 Pull the end of the paper past the recorder platen.**
- 5 While holding the recorder door open, (to the right) press the platen down over the paper.**
- 6 Allow the door to close over the platen roller.**

Cleaning the Recorder Printhead

If the ECG strip has light or varying density printing, clean the printhead to remove any possible buildup of paper residue.

To clean the printhead, perform the following steps.

- 1 Slide the recorder door to the right of the defibrillator. The paper platen will tilt up.**
- 2 Remove the paper roll.**
- 3 Clean the printhead surface (above the brush) with rubbing alcohol and a cotton swab.**
- 4 Re-install the paper roll.**
- 5 While holding the recorder door open, (to the right) press the platen down over the paper.**
- 6 Allow the door to close over the platen roller.**

Maintaining the Battery

The sealed lead-acid battery used in the CodeMaster XL+ will provide optimum life when the unit is continuously connected to AC power and fully charged after each use. To fully charge a depleted battery requires 18 hours of continuous charge time. Because it is not always possible to allow a full charge cycle between uses, the CodeMaster XL+ was designed to charge a depleted battery to 90% of its capacity in approximately two hours. However, battery capacity and battery life will be reduced if the battery is not allowed to fully charge after each use. For improved battery life, applications where the CodeMaster XL+ is used frequently between full charge cycles, consider ways to decrease the number of instrument uses between full charge cycles.

When the instrument is not plugged into AC power, some current is drawn from the battery to maintain memory and startup logic. If the battery is to be stored for extended periods (more than one month) without AC power, observe the caution in Appendix A, *Installing and Charging the Battery*.

This battery was selected because it provides optimum performance and battery life over a wide range of operating conditions. The life expectancy of this battery is dependent on many variables, including temperature and usage. Periodically check the battery capacity to determine whether to replace it. The battery capacity check is described in the next section.

NOTE



When plugged into AC power the CodeMaster XL+ will function normally with no battery installed, however the time required to charge the defibrillator will increase.

WARNING



If the CodeMaster XL+ is operated without a battery installed, clearly mark that the instrument does not have a battery and requires AC power to operate. When a CodeMaster XL+ has no battery installed and is plugged into AC power, the front panel **AC POWER** light will be on, and the **BATT CHRG** light will be off.

Battery Capacity Check

Perform the battery capacity check at least once every six months. If the battery is frequently depleted without adequate time for full charge cycles, perform the check more often.

A new battery will provide a minimum of 2.5 hours of monitoring time. Hewlett-Packard recommends that you replace the battery when it fails to provide a minimum of 2.5 hours of monitoring time or 10 minutes of Low Battery warning time when starting from a fully charged battery.

To test the battery capacity, perform the following steps. This test will require up to 3.5 hours to perform (not including battery charge times). Allow 30 hours to test the battery and have it ready for use again.

- 1 Charge the battery by plugging the instrument into AC power for eight hours. Verify that the **AC POWER** and **BATT CHRG** lights are on.
- 2 Verify that there is recorder paper in the recorder.
- 3 Depress **Sync** **HR Alarm** while turning the Energy Select control from Off (Standby) to Monitor On. While pressing **Sync**, turn the Energy Select control from the Off position to the Monitor On position. The Setup/Diagnostic Service menu will appear in a moment.
- 4 Press **▼ ECG Size** to highlight the Test Battery diagnostic test.
- 5 Unplug the instrument from AC power (the **AC POWER** and **BATT CHRG** lights are off).
- 6 Press **Lead Select** to start the Battery Capacity Test.
 - The starting battery voltage, the current battery voltage, the Monitor elapsed time, and the Low Battery warning elapsed time will be displayed.
 - When the test is finished, the recorder will print the final values of the displayed results just prior to turning off the instrument.
 - Replace the battery if the recorded value for elapsed monitor time is less than 2.5 hours or the value for elapsed Low Battery warning time is less than 10 minutes.

- 7 Turn the Energy Select control to the Off (Standby) position.
- 8 Fully recharge the battery by plugging the instrument into AC power for 18 hours. Verify that the **AC POWER** and **BATT CHRG** lights are on.
- 9 The instrument is now ready to be returned to service.

CAUTION



The battery can be permanently damaged if left uncharged for a prolonged period.

Replacing the Battery

To install a new battery, perform the following steps:

- 1 Remove AC power from the instrument.
- 2 Place the instrument upside down on a workbench.
- 3 Turn the two battery compartment quarter-turn screws counter-clockwise.
- 4 Lift out the battery door.
- 5 Disconnect the battery connector and remove the battery.
- 6 Connect the new battery and slide it into the compartment.
- 7 Replace the battery door and secure it by turning the two retaining screws a quarter-turn clockwise.
- 8 Perform the Battery Capacity Check above before placing instrument in service.

Cleaning Exterior Surfaces

The CodeMaster XL+ and its accessories are chemically resistant to common hospital cleaning solutions and non-caustic detergents. The following list includes some approved cleaning solutions.

- 90% Isopropyl alcohol (except adaptors and patient cable)
- Soap and water
- Chlorine bleach (30 ml/l water)
- Ammonia-based cleaners
- Keep the outside of the instrument clean and free of dust and dirt. Clean the paddles thoroughly to prevent build-up of dried electrolyte paste.
- Do Not allow any fluids to penetrate the instrument case. Avoid pouring fluid on the unit while cleaning.
- Do Not use abrasive cleaners, or strong solvents such as acetone, or acetone-based compounds.
- Clean the display screen carefully. It is especially sensitive to rough handling and subject to scratching.
- Do not steam sterilize the monitoring leads, submerge them for prolonged periods, or heat them above 50°C. If metallic surfaces become oxidized, clean them with a very light abrasive (toothpaste). Do not use highly abrasive cleaners such as steel wool or silver polish.
- Do not steam or gas sterilize the external paddle set.

Cleaning and Sterilizing the Internal Paddles

- Clean the electrode surface and handle with standard hospital solution.
- Do not use acetone and ammonia-based cleaners.
- Use a small, soft brush with cleaning solution to clean any contamination from the electrode surface and edges.
- Before sterilizing, clean any excessive residue which accumulates on the handles or electrode surfaces.

Internal paddle sets can be sterilized using any of four different methods. Three of these methods use steam sterilization, the fourth method uses Ethylene Oxide (EtO). Although the internal paddle sets are constructed using the highest quality materials, the severe conditions of steam sterilization will limit their useful life.

Hewlett-Packard Company provides the standard, 90 day *consumables* warranty on the paddles. Since the useful life of these paddles is limited, replace them when functionality or appearance is questionable. When using the sterilization procedures described in this guide, the paddles will withstand approximately 200 sterilization cycles. This will vary depending on the equipment and the process used.

Steam Sterilizing the Internal Paddles

Prevacuum Sterilization

Preparation	Double wrapped in sterilization grade wrap
Exposure Temperature	132–135°C (270–275°F)
Exposure Time	3 minutes

Flash Sterilization

Preparation	Unwrapped
Exposure Temperature	132–135°C (270–275°F)
Exposure Time	3 minutes

Gravity Sterilization

Preparation	Double wrapped in sterilization grade wrap
Exposure Temperature	121–123°C (250–254°F)
Exposure Time	30 minutes

Ethylene Oxide Sterilization

Temperature	54°C ± 1°C (130°F ± 2°F)
Relative Humidity	60% ± 20%
EtO Concentration	600 mg/L ± 30 mg/L 12/88 EtO
Exposure Time	1 hour 45 minutes
Aeration Time	18 hours

Supplies

The Hewlett-Packard warranty is only assured if you use Hewlett-Packard approved accessories and replacement parts.

Battery

M1758A	Battery assembly
---------------	------------------

Patient Cables

M1731A	Patient cable 6-pin/3-wire AHA
M1732A	Patient cable 5-wire/6-pin AHA
M1733A	Patient cable 8-pin/3-wire AHA
M1734A	Patient cable 5-wire/8-pin AHA
M1735A	Patient cable 8-pin/3-wire IEC
M1736A	Patient cable 5-wire/8-pin IEC
M1500A	Trunk cable 3-wire/12-pin - AHA
M1605A	Lead set 3-wire - AHA
M1510A	Trunk cable 3-wire/12-pin - IEC
M1615A	Lead set 3-wire - IEC
M1520A	Trunk cable 5-wire/12-pin - AHA
M1625A	Lead set 5-wire - AHA
M1530A	Trunk cable 5-wire/12-pin - IEC
M1635A	Lead set 5-wire - IEC

Power Cords

8120-1692	Continental Europe
8120-4464	Australia, New Zealand
8120-1703	United Kingdom, Ireland
8120-2957	Denmark
8120-1692	Italy
8120-5213	North America
8120-2296	Switzerland
8120-4600	South Africa

Paddles

M1740A	Internal paddle adaptor
M1740B	Internal paddle adaptor (yellow)
M1741A	Switchless Internal paddle set 7.5cm diameter
M1742A	Switchless Internal paddle set 6.0cm diameter

**Maintaining the Defibrillator
Supplies**

M1743A	Switchless Internal paddle set 4.5cm diameter
M1744A	Switchless Internal paddle set 2.8cm diameter *
M1784A	Switched Internal paddle set, 7.5 cm diameter
M1785A	Switched Internal paddle set, 6.0 cm diameter
M1786A	Switched Internal paddle set, 4.5 cm diameter
M1787A	Switched Internal paddle set, 2.8 cm diameter *
M1746A	Anterior/anterior paddle set - PCI
M1746B	Anterior/anterior paddle set (yellow) - PCI
M1747A	Anterior/anterior paddle set - no PCI
M1747B	Anterior/anterior paddle set (yellow) - no PCI
M1748A	Anterior/anterior paddle set, EtO sterilizable - no PCI

*Does not comply with IEC 601-2-4 standard for contact area.

Pads

M1749A	Defibrillation/pacing pads
M1750A	Pads adaptor
M1750B	Pads adaptor (yellow)
M1781A	Pads adaptor test load

Electrolyte (Redux®)

651-1008-050 Redux®

Paper

40457C/D Recorder paper

Cases

M1778A	Carry case
M1779A	Accessory pouch

**Maintaining the Defibrillator
Supplies**

Wall Mount Kit

M1722-80001 **Wall mount hardware**

**Maintaining the Defibrillator
Supplies**

Installation and Setup

Installation

The CodeMaster XL+ is ready for operation when the following tasks have been properly performed:

- Install battery.
- Charge battery (for 24 hours).
- Install paper.
- Make sure that the paddle set connector is seated and locked.
- Select configuration settings; set date and time.

Line Voltage Settings

The defibrillator automatically adjusts to the line voltage that is supplied (from 100–230 VAC \pm 15% at 50/60 Hz). No manual setting or adjustment is required.

WARNING



Use only three wire power cords with three-pronged grounded plugs. Make sure that the outlet accepts the three-pronged plugs and is grounded. Never adapt a three-pronged plug to fit a two-pronged outlet.

Installing and Charging the Battery

To install the battery, refer to the battery replacement procedure in “Replacing the Battery” (page 8-5). After installing the battery, connect the power cord to the back of the defibrillator, then plug the cord into an AC outlet. The green **AC POWER** and **BATT CHRG** indicators on the front panel should light up. (The **AC POWER** indicator lights when the instrument is plugged into AC power; the **BATT CHRG** indicator is on when the battery is installed and the instrument is plugged into AC power.)

NOTE



To ensure full battery capacity, charge the battery for 24 hours following its installation in the defibrillator.

The defibrillator operates from either battery or AC power.

Use only HP battery assembly M1758A.

WARNING




To avoid the possibility of hazardous electrical shock, unplug the instrument from the AC power source before installing or replacing the battery.

CAUTION



If the defibrillator will be stored for longer than one month without AC power, first charge the battery for 48 hours, then remove the battery from the unit. Note on the instrument that the battery has been removed. Store the battery in a cool, dry location. Recharge a stored battery for at least 24 hours every six months. This will ensure that the battery does not completely discharge while in storage. The battery's shelf life is longer with cooler temperatures, but the battery must not be stored below freezing. After an extended storage period, the battery should be tested using the "Battery Capacity Check" (page 8-4).

Figure 1-1



Installing the Battery

Loading the Recorder Paper

The defibrillator recorder uses two-inch wide, thermal paper (HP 40457C/D). To load the paper, refer to the procedure in “Changing the Recorder Paper” (page 7-1).

Connecting Paddles and Patient Cables

The defibrillator has a paddles connector for attaching pads/paddles sets and an ECG Input connector for attaching leads.


Paddles Connector

The defibrillator connector accepts external paddles, external adhesive pads, or internal paddles.

Connecting Paddles or Pads To connect external paddles, internal paddles, or adhesive pads to the defibrillator, perform the following steps:

- 1** Slide the paddle connector lock on the paddles plug to the **unlock** position. To do this, push the lock towards the top of the connector.
- 2** Insert the paddles/pads adaptor cable plug into the paddles connector on the defibrillator, as shown in Figure 7-1 on page 7-1.
- 3** Slide the paddle connector lock to the **lock** position, to latch the plug in place.

Figure 1-2



Connecting External Paddles, Adhesive Pads, or Internal Paddles

ECG Input Connector

The ECG Input connector on the defibrillator is a 6-, 8-, or 12-pin connector, depending on the option purchased with the instrument. For each connector option, several different patient cables can be used for various ECG sources and applications.

Refer to Table 3-1, "Cardiac Monitoring Configurations," on page 3-1 for a list of available patient cables and lead sets, and their part numbers.

NOTE



3 wire = RA, LA, LL
5 wire = RA, RL, LA, LL + V (C)

Connecting a Patient Cable The 3-wire or 5-wire patient cable connects to the ECG Input connector located on the front of the defibrillator, behind the carrying handle. The patient cable plug has 6-, 8-, or 12-pins. Before installing the patient cable, make sure that the pin count of the patient cable plug matches the pin count of the ECG Input connector. To install the patient cable:


- 1 Align the keyed cable plug with the slot in the ECG Input connector. See Figure 1-3 on page A-7.
- 2 Push the cable plug firmly into the ECG Input connector.

NOTE



When a lead is selected for monitoring, the message **LEADS OFF** appears on the display if the patient cable falls off or is incorrectly connected. Also, a dashed line appears on the display in place of an ECG trace.

Figure 1-3



Connecting a Patient Cable

Setup

To configure the CodeMaster XL+ you must use the setup menus. Table 1-1, "Setup Menu 1 Settings," on page A-9 and Table 1-2, "Setup Menu 2 Settings," on page A-10 list the choices you can make on the setup menus. Perform the following steps to configure the CodeMaster XL+.

- 1 Depress **Sync** **HR Alarm** while turning the Energy Select control from Off (Standby) to Monitor On. A menu will appear with the following choices.
CALIBRATE DEFIB
SETUP MENU 1
SETUP MENU 2
RESTORE FACTORY SETTINGS
PRINT LOG
TEST DEFIB
TEST ECG
TEST CRT
TEST RECORDER
TEST CONTROLS
TEST INDICATORS
TEST BATTERY
TEST PACER
- 2 Select SETUP MENU 1 by pressing **ECG Size** until the highlight appears on **SETUP MENU 1**.
- 3 Press **Lead Select**. SETUP MENU 1 will appear with the current setup values displayed.
- 4 Press **ECG Size** or **ECG Size** until the highlight appears on the value you wish to change.
- 5 Press **Lead Select**.
- 6 Press **ECG Size** or **ECG Size** to scroll through the available choices for this parameter.
- 7 When the choice you want is displayed, press **Lead Select** to set your choice.
- 8 Press **ECG Size** or **ECG Size** until the highlight appears on the next value you wish to change.
- 9 Repeat steps 5 through 8 until you are finished configuring the settings in setup menu 1.

10 Depress both sides of the  key at once to return to the setup/diagnostic menu.

NOTE



You must depress both sides of the  key at the same time to return to the setup/diagnostic menu.

11 To change settings in setup menu 2, select **SETUP MENU 2** from the setup/diagnostic menu and repeat the above steps.

12 Turn the Energy Select control to Off (Standby) to leave setup/diagnostic mode.

You can use the factory setting for most values by selecting **RESTORE FACTORY SETTINGS** from the setup/diagnostic menu.

WARNING



The setting values can have critical impact on how your defibrillator operates.

If the CodeMaster XL+ loses your configuration settings, it will display the message **SETUP LOST** on the screen and use the factory settings for all setup values. To clear the **SETUP LOST** message, go to setup menu 1 or setup menu 2. The "Setup Lost" condition will be cleared when you view the setup menus. You do not have to change the values if you want to keep the factory settings.

Table 1-1, "Setup Menu 1 Settings," on page A-9 and Table 1-2, "Setup Menu 2 Settings," on page A-10 show configurable parameters on the instrument.

Table 1-1 Setup Menu 1 Settings

Setting	Choices	Description
Language	English, Dutch, Swedish, French, German, Italian, Spanish, Finnish, Danish, Norwegian	Printed and displayed text language
Upper Alarm Limits (UAL)	120, 140, 160, 20 to 280 in increments of 5	Upper heart rate limits LAL to 280

Table 1-1 Setup Menu 1 Settings

Setting	Choices	Description
Lower Alarm Limits (LAL)	40, 60, 90, 20 to 280 in increments of 5	Lower heart rate limits 20 to UAL
Time	HH:MM	Current time
Date	DD MMM YY	Current date
Armed Tone	ON, OFF	Beep on Charge done
CRT Alerts	ON, OFF	Beep on alert message
Alert Volume	3 to 15	Alerts volume
Mode after CV	SYNC, DEFIB	Specifies mode after cardioversion
Pacer Rate (pacer option only)	70 (ppm), 40 to 180 in increments of 10	Sets the pacer rate and the initial power-on default rate.
Pacer Output (pacer option only)	30 (mA), 30 to 200 in increments of 10	Sets the pacer current and the initial power-on default pacer current.

Table 1-2 Setup Menu 2 Settings

Setting	Choices	Description
Recorder Delay	Delay 6s, No Delay	Printout is delayed 6 seconds or immediate
Recorder BW	Monitor, Diagnostic	Bandwidth (recorder only)
Record on Mark	ON, OFF	Records during mark
Record on Charge	ON, OFF	Records during charge
Record on Shock	ON, OFF	Records during discharge
Record on Alarms	ON, OFF	Records during alarms.
Post Shock Data	ON, OFF	ON prints the delivered energy statistics. OFF prints Energy Select control setting.
Power On Lead	PADDLES, LEAD I, LEAD II, LEAD III	Sets the ECG monitoring source that appears when you turn on the instrument.
Patient Cable	3 WIRE, 5 WIRE	

Table 1-2 **Setup Menu 2 Settings**

Setting	Choices	Description
Notch Filter	60 HZ, 50 HZ, ON, OFF	
ECG Trace	SWEEP, SCROLL	ECG trace style

Specifications

Defibrillator

Waveform: Damped sinusoidal (Lown).

Output Energy (Delivered): 2, 3, 5, 7, 10, 20, 30, 50, 70, 100, 150, 200, 300, and 360 joules.

Charge Control: Push-button on apex paddle and on front panel.

Charge Time: Less than 5 seconds to 360 joules with battery present. Less than 15 seconds to 360 joules on AC only.

Armed Indicators: Charge done tone, charge done lamp on apex paddle, and available energy indicated on display.

Paddle Contact Indicator: 3-color LED bar graph array on STERNUM paddle indicates quality of defibrillator paddle contact before discharge.

Paddles: Standard paddles are anterior/anterior, adult and pediatric. Adult electrodes (83 cm sq) slide off to expose pediatric electrodes (21 cm sq). Paddle cord is 10 ft (3 m). Full range of internal paddles available.

Synchronizer: SYNC message appears on monitor and is annotated periodically on recorder while in synchronous mode. An audible beep sounds with each detected R-wave, while a marker on the monitor and sync designator on the recorder strip indicate the discharge point.

Environmental Operating Conditions: 0 to 55 deg C, 15 to 95% relative humidity, 15,000 ft altitude.

Environmental Storage Conditions: -20 to 70 deg C, 90% relative humidity for 24 hours at 65 deg C, 15,000 ft altitude.

Monitor

Inputs: ECG may be viewed through paddles or patient cable. Lead I, II, III, or PADDLES selectable. Additional leads (avR, avF, AvL, V Leads) and PADS are available. Monitor and recorder indicate selected ECG source.

Lead Fault: LEADS OFF message and dashed baseline appear on monitor if a lead becomes disconnected.

Common Mode Rejection: Greater than 100 dB measured as per AAMI standards for cardiac monitors (EC13).

Display Size and Type: 5 inch (12.7 cm) diagonal CRT for 4 seconds of ECG data on screen; non-fade, fixed trace. Scrolling trace is selectable.

Sweep Speed: 25 mm/sec nominal.

Frequency Response: 0.5 to 40 Hz.

Heart Rate Display: Digital readout on monitor from 15 to 300 BPM.

Heart Rate Alarms: Three configurable pairs of high and low heart rate alarm limits from 20 to 280 BPM.

ECG Output: 1V/mV.

Patient Cable Length: 10 ft (3 m).

Thermal Array Recorder

Event Summary: Stores and prints 3 seconds pre- and 8 seconds post-critical event data for up to 28 events. Data retained after unit is turned off.

**Installation and Setup
Specifications**

Annotates: Time, date, HR, event marker, ECG mode, defibrillator mode, selected energy, actual delivered energy, peak current, and patient impedance.

Speed: 25 mm/sec.

Paper Size: 50 mm by 30 m (100 ft).

Recorder Mode: Automatically documents events and ECG during defibrillation episodes. The recorder can be configured to run in either real time or with a six-second delay.

Frequency Response: 0.5 to 40 Hz or 0.05 to 150 Hz selectable.

Size and Weight

Dimensions: 7.9 in H by 11.8 in W by 15.6 in L (20 cm H by 30 cm W by 39.7 cm L).

Weight: 24 lbs (10.9 kg). Includes external paddles, battery, and recorder paper.

Battery

Type: Rechargeable sealed lead-acid battery. 4 Ah, 12 V nominal.

Charge Time: 2 hours to 90% of full capacity. 18 hours to 100% capacity. Repeated charging to less than 100% will reduce useful life of battery.

Capacity: 2 1/2 hours monitoring or fifty (50) full-energy discharges or 1 hour monitoring and recording.

Battery Indicators: Illuminated LED indicates battery is charging. LOW BATTERY message appears on monitor when limited battery capacity remains.

AC Input: 100 to 230 VAC +/-15%, 50 to 60 Hz.

External Pacemaker (Optional)

Current Pulse Amplitude: 30 mA to 200 mA.

Pulse Width: 20 msec.

Rate: 40 ppm to 180 ppm.

Modes: Demand or fixed rate.

Refractory Period: 40 to 80 ppm 340 msec; 90 to 180 ppm 240 msec.

SpO₂ (Optional)

Measurement Range: 0 to 100%.

Accuracy with HP M1190A sensor: 1 standard deviation, 65-80%: $\pm 2.5\%$, 80-100%: $\pm 1.5\%$, resolution: 1%.

Averaging: 8 beats.

SpO₂ alarm limits — range: 100/90, 100/85, 100/80.

SpO₂ alarm delay: ten seconds after value drops below the low alarm setting.

INOP alerts: Triggered by disconnected sensor, noisy signal, light interference or low signal.

Pulse rate measurement range: 30 to 300 bpm $\pm 1\%$; resolution 1 bpm.

Pulse amplitude indicator: Indicates pulsatile activity.

Calling for Service

United States

Eastern Region
Tel: 1-800-245-4406

Western Region
Tel: 1-800-428-0949

Southern Region
Tel: 1-800-845-8374

Midwestern Region
Tel: 1-800-323-2280

Canada

Eastern Region
Tel: 1-800-361-9790

Central Region
Tel: 1-800-387-3900

Western Region
Tel: 1-800-661-5626

Other International Areas

European Headquarters
Hewlett-Packard S.A.
150 Route du Nante d'Avril
P.O. Box CH-1217, Meyrin 2
Geneva, Switzerland
Tel: 22 838-111

France
Hewlett-Packard France S.A.
Parc d'activite du Bois Briard
2, avenue du Lac
F-91040 Evry Cedex
Tel: (1)60 77 83 83

**Austria, Eastern Europe, Yugosla-
via**
Hewlett-Packard Ges.m.b.H
Lieblgasse 1
P.O. Box 72
A-1222 Vienna
Tel: 222/25 00-0

Germany
Hewlett-Packard GmbH
Hewlett-Packard-Strasse
D-6380 Bad Homburg
Tel: 0 61 72/400-0

Belgium
Hewlett-Packard Belgium SA/NV
Boulevard de la Woluwe, 100
Woluwedal
B-1200 Brussels
Tel: 2/761 31 11

Mediterranean and Middle East
Hewlett-Packard S.A.
Atrina Center
32, Kifissias Ave
Paradissos-Amaroussion GR-Athens
15125
Tel: 1/682 88 11

Other International Areas

Netherlands

Hewlett-Packard Nederland B.V.
Startbaan 16
1187 XR Amstelveen
P.O. Box 667
NL-11 80 AR Amstelveen
Tel: 20/547 69 11

Norway

Hewlett-Packard Norge A/S
P.O. Box 34
Osterdalen 16-18
N-1345 Osteras
Tel: 2/24 60 90

Denmark

Hewlett-Packard A/S
Kongevejen 25
DK-3460 Birkerød
Tel: 2/81 66 40

Finland

Hewlett-Packard OY
Piispankalliontie 17
02200 Espoo
Tel: 0/887 21

Spain

Hewlett-Packard Espanolas S.A.
Ctra. de la Coruna
km 16, 400
Las Rozas
E-Madrid
Tel: 1/637 00 11

Sweden

Hewlett-Packard Sverige AB
Skalholtskatan 9, Kista
Box 19
S-163 93 Spanga
Tel: 8/750 20 00

United Kingdom

Hewlett-Packard Limited
Harman House
No. 1 George Street
Uxbridge Middlesex UB8 1YH
Tel: 895/720 20

Italy

Hewlett-Packard Italiana S.p.A.
Via G. di Vittorio, 9
20063 Cernusco S/N (MI)
Tel: 2/92 36 91

Index

Numerics

- 3-wire patient cable, 3-1, 3-2, 8-9, A-6
- 50J MAXIMUM, 2-7, 2-8, 7-2
- 5-wire patient cable, 3-1, 3-2, 3-3, 8-9, A-6

A

- AC POWER, 2-4, 4-5, 7-7, 7-9
- AC power, 1-7, 1-12, 4-5, 7-9, 8-5
 - operation, 1-10
 - pacing, 5-2
- adaptor cable
 - internal paddles, 2-6
- adhesive pads, 2-5, 5-1
- adjusting
 - ECG size, 3-7, 4-3
 - heart rate alarm limits, 3-7
 - QRS beeper volume, 3-7
- alarm
 - battery, 1-11
 - heart rate, 1-4
- alarms violation
 - recordings, 3-10
- alcohol
 - for cleaning, 8-6
- Alert Volume, A-10
- Armed Tone, A-10
- arrhythmias
 - synchronized cardioversion as treatment, 4-1
- artifact, 4-2
- asynchronous events, 3-9
- ATTACH PADS, 5-3
- audible indicators, 1-8
- AUTOGAIN, 7-7
- Autogain, 3-7, 4-3
- automatic recordings, 3-9
- available energy, 2-3, 3-10

B

- BATT CHRГ, 1-7, 1-10, 2-4, 4-5, 7-7, 7-9
- battery
 - charging, A-1
 - charging the, 1-11
 - compartment, A-3
 - connector, A-3
 - discharges, 1-11
 - lead-acid, 8-3
 - life, 1-11
 - operates without a, 1-12
 - plug, A-3
 - replacing, 8-5, A-1

- battery alarm, 1-11
- battery capacity, 1-11
 - check, 8-4
 - of new, 8-4
- battery door, 8-5
- battery operation, 1-10
- bell symbol, 3-7

C

- cable
 - 3-Wire, 3-1
 - 3-wire, 3-2, 8-9
 - 5-wire, 3-1, 3-3, 8-9
 - pads adaptor, 2-5, 4-2
 - patient, 3-1
- calibration pulse, 3-9
- capturing a beat, 5-4
- cardioversion, 3-5
 - elective, 3-1
 - performing, 4-2
 - synchronized, 3-1
- caution, iv
- changing the paper, 8-1
- charge, 2-2
- Charge button, 1-2, 4-3, 5-4
- Charge Done
 - indicators, 2-3, 4-5, 5-4, 7-10
 - light, 1-7
 - tone, 1-8
- charge done
 - indicators, 4-3
- charging
 - battery, A-1
 - charging the battery, 1-11
- CHECK RECORDER, 3-10, 7-3
- checking
 - battery capacity, 8-4
 - checking connections, 3-6
 - checks
 - battery capacity, 8-4
 - operational, daily, 7-9
 - operational, every shift, 7-9
 - chlorine bleach
 - for cleaning, 8-6
 - cleaning
 - exterior, 8-6
 - paddles, 8-6
 - recorder printhead, 8-2
 - cleaning solutions, 8-6
 - clearing the Event Summary, 3-8
 - CodeMaster XL+ defibrillator/monitor, 1-1
 - configuring

- set date and time, A-1
- settings, A-1
- configuring the defibrillator/monitor, A-7
- connecting
 - patient cable, A-7
- connections
 - checking, 3-6
 - pacing pad, 5-3
- connector
 - defibrillator, A-3
 - lock, 2-5, A-5, A-7
 - paddle, 2-5
- controls
 - defibrillator, 1-2, 2-1
 - monitor, 1-4
 - pacer, 1-5
 - recorder, 1-3
 - SpO₂, 1-6
- conventions, manual, iv
- countershock, 4-3, 4-5
- CRT Alerts, A-10
- CRT alerts, 1-8

D

- Date, A-10
- decreasing energy level, 2-3
- DEFIB DISARMED, 7-3
- DEFIB FAILURE, 7-1
- defibrillating
 - adhesive pads, through, 2-5
 - during pacing, 5-4
 - external pads, through, 2-5
 - internally, 2-6, 2-8
 - three steps, 2-1
 - through pacer pads, 5-1
 - with alternate paddles, 2-4
 - with pediatric paddles, 2-4
- defibrillator
 - after using, 2-4, 4-5
 - charging, 2-2
 - connector, A-3
 - controls, 1-2, 2-1
 - disarming, 1-8
 - discharging, 2-3
 - if no shock, 4-5
 - maintaining the, 8-1
 - operation check, 4-3
 - removing power from, 1-8
 - replacing the battery, 8-5
 - turning on the, 2-2
- defibrillator connector
 - external adhesive pads, A-4

Index

- external paddles, A-4
 - internal paddles, A-4
 - uses, A-4
 - defibrillator mode, 1-7, 4-5
 - defibrillator protection, iv
 - delayed recorder mode, 3-10
 - Delivered Energy, 2-3
 - delivered energy, 4-3
 - delivered energy check, 7-9
 - Delivered Energy Display, 5-4
 - demand mode, 1-5, 5-1, 5-2
 - diagnostic mode, A-9
 - diagnostics, 7-7
 - diastole
 - pacemaker pulse position in diastole, 5-3
 - disarm
 - recordings, 3-10
 - disarming
 - automatic, 1-9
 - charged defibrillator, 1-8
 - discharge
 - stored energy, 4-4
 - E**
 - ECG gain mode, 7-7
 - ECG Input
 - connecting, A-7
 - ECG Input connector
 - options, A-5
 - uses, A-5
 - ECG output jack, 4-1
 - ECG Size, 3-9, 4-3
 - ECG size
 - adjusting, 3-7, 4-3
 - ECG Size button, 1-4, 3-7
 - ECG source, 1-4, 4-2
 - selecting, 3-6
 - ECG strip, 3-8, 4-3
 - annotating, 1-3
 - light printing, 8-2
 - printing, 1-3
 - reading the, 3-9
 - elective cardioversion, 3-1
 - electrode, 3-1, 7-9
 - monitoring, 5-2
 - placing 3-wire, 3-2
 - placing 5-wire, 3-3
 - electrolyte paste, 7-9
 - applying, 2-2
 - energy
 - available, 2-3
 - decreasing, 2-3
 - delivered, 2-3
 - increasing, 2-3
 - resetting, 2-3, 4-4
 - selecting, 1-2, 2-2
 - Energy Select control, 1-2, 1-8, 2-2, 4-5
 - ES, 1-4, 3-8
 - Event Summary, 1-4, 2-4
 - clearing, 3-8
 - printing, 3-8
 - reviewing, 3-8
 - stopping, 3-8
 - events, 3-8
 - external adhesive pads
 - defibrillator connector, A-4
 - external monitor
 - using with the, 4-1
 - external paddles, 2-4
 - connecting, A-4
 - defibrillator connector, A-4
 - external paddles adaptor cable
 - connecting, A-4
 - connecting to defibrillator, A-4
 - external pads, 2-4, 4-2
 - adaptor, 2-5
 - defibrillating through, 2-5
 - monitoring through, 2-5
 - F**
 - factory settings, A-9
 - fixed mode, 1-5, 5-1, 5-2
 - forming leads, 3-3
 - functional inspection, 7-7
 - G**
 - gain bar, 3-7
 - H**
 - Heart Rate, 5-4
 - heart rate alarm, 1-4, 1-8
 - adjusting limits, 3-7
 - pacemaker patients, 1-10
 - reviewing limits, 3-7
 - symbol, 3-7
 - violated limits, 3-7
 - heart rate alarms
 - during pacing, 5-4
 - Heart Rate Alarms violation, 3-8
 - high voltage
 - paddles, 1-9
 - HR Alarm button, 1-4, 3-7
 - I**
 - increasing energy level, 2-3
 - indicator
 - paddle contact, 4-5
 - indicator lights, 1-7
 - indicators
 - audible, 1-8
 - Charge Done, 2-3, 4-5, 5-4, 7-10
 - inspection
 - functional, 7-7
 - installation
 - required tasks, A-1
 - installing
 - battery, A-1
 - paddle set, A-1
 - paper, A-1
 - the defibrillator, A-1
 - inter paddles
 - switched, 2-8
 - internal, A-4
 - internal defibrillation, 2-6, 2-8
 - energy limit, 2-8
 - energy limits, 2-7
 - internal paddle adaptor, 8-9
 - internal paddle set, 8-9
 - internal paddles, 2-4, 2-6
 - internal paddles adaptor, 2-6
 - internal paddles adaptor cable, 2-6
 - internal paddles, switched, 2-6
 - internal paddles, switchless, 2-6
- L**
 - Language, A-9
 - Lead Select button, 1-4, 4-2, 5-3, A-8
 - leads
 - forming, 3-3
 - preparing for monitoring, 3-2
 - LEADS OFF, 3-6, 5-3, 5-4, 7-2, 7-7
 - lights, 1-7
 - AC POWER, 4-5
 - BATT CHRG, 4-5
 - SYNC, 7-3
 - line voltage
 - automatic setting, A-1
 - loading recorder paper, A-3
 - locking
 - paddle set connector, A-1
 - LOW BATTERY, 1-11, 7-2
 - Low Battery warning time, 8-4
 - Lower Alarm Limits, A-10
- M**
 - maintaining
 - defibrillator, 8-1
 - Mark button, 1-3, 3-8, 7-7
 - marker dot

Index

- on R-wave, 4-3
- on T-wave, 4-3
- menu
 - SETUP MENU 1, A-9
 - setup/diagnostic, A-9
- mode
 - Autogain, 4-3
 - defibrillator, 1-7, 4-5
 - delayed recorder, 3-10
 - Demand pacer, 5-2
 - demand, pacing, 1-5, 5-1
 - diagnostic, A-9
 - ECG gain, 7-7
 - Fixed pacer, 5-2
 - fixed, pacing, 1-5, 5-1
 - non-delayed recorder, 3-10
 - paddles ECG, 7-7
 - setup, A-9
 - sync, 1-2, 1-3, 4-5
 - synchronized shock, 1-7
 - synchronous, 3-5
- Mode button, 1-5, 5-3
- MONITOR FAILURE, 7-1
- Monitor Operating Controls, 1-4
- monitoring, 1-4
 - adhesive pads, through, 2-5
 - cables, 3-1
 - leads, 3-1
 - on AC power, 3-1
 - on battery, 3-1
 - preparing the leads, 3-2
 - with paddles, 3-5
- monitoring electrodes
 - separate from pacing pads, 5-2
- monitoring equipment, 4-1
- N**
- NO PADDLES, 5-2, 7-2
- NO PADS, 7-5
- non-delayed recorder mode, 3-10
- O**
- operational checks, 7-8
- Output button, 1-5, 5-4
- output current
 - pacing, 1-5
- oxygen saturation, 6-1
- P**
- pacer
 - control panel, 5-1
 - controls, 1-5
 - messages, 7-5
 - output, 5-2
- PACER FAILURE, 7-5
- Pacer On, 1-5
- Pacer On button, 5-2
- Pacer Output, A-10
- pacer parameters, 5-2
- pacer pulse
 - delivery to patient, 5-3
 - positioning in diastole, 5-3
- Pacer Rate, A-10
- pacer rate, 5-2
- PACER STOP, 5-2, 5-4
- PACING, 5-3
- pacing, 3-1
 - defibrillating during, 5-4
 - demand mode, 1-5
 - external transcutaneous, 5-1
 - fixed mode, 1-5
 - on AC power, 5-2
 - on battery, 5-2
 - output current, 1-5
 - with electrodes, 5-1
- pacing pad, 5-4
- paddle connector, 2-5
- paddle connector lock, 2-5
- Paddle Contact Indicator, 1-7
- paddle contact indicator, 4-5
- paddle placement
 - cardioversion, 4-4
- paddle set, 8-9
- paddle, Sternum
 - placement, 2-2
- paddle, Apex
 - placement, 2-2
- paddles
 - alternate, 2-4
 - applying the, 2-2
 - external, 2-4
 - installing, A-3
 - internal, 2-4, 2-6
 - internal, switched, 2-6, 2-8
 - internal, switchless, 2-6
 - monitoring with, 3-5
 - pediatric, 2-4
 - plug, A-5
 - preparing the, 2-2
 - pressure, 2-2
 - removing from holders, 2-2
- paddles connector
 - uses, A-4
- paddles/pads
 - installing, A-4
- paddles/pads connector
 - locking, A-4
- pads
 - external, 4-2
 - external adhesive, 5-1
 - multifunction, 5-5
 - pacing, 5-2, 5-4
- pads adaptor cable, 2-5, 3-2, 4-2
- PADS OFF, 2-5, 3-6, 4-2, 5-2, 5-4
- paper
 - loading, A-3
- paper platen, 8-1, 8-2
- patient cable, 8-9
 - 3-wire, 3-2
 - 5-wire, 3-2, 3-3
 - connecting, A-6, A-7
 - keyed plug, A-7
- patient cables
 - installing, A-3
- pediatric defibrillation, 2-4
- pediatric paddles, 2-4
- placing electrodes
 - 3-wire, 3-2
 - 5-wire, 3-3
- placing paddles for cardioversion, 4-4
- Post CV Mode, A-10
- Post Shock Data
 - setting, 3-10
- power cord
 - worn, 7-9
- preparing the paddles, 2-2
- printhead
 - cleaning, 8-2
- Printing
 - Event Summary, 3-8
- printing
 - calibration pulse, 3-9
 - current ECG status, 3-9
 - ECG strip, 1-3
 - monitor status, 3-9
 - shock delivery statistics, 3-10
- pulse oximetry, 6-1
- Q**
- QRS beeper, 1-8
- QRS beeper volume, 3-7
- R**
- Rate, 1-5
- Rate button, 5-3
- Record button, 1-3, 3-9, 7-7
- recorder
 - changing the paper, 8-1
 - cleaning the printhead, 8-2

Index

- loading paper, A-3
- recorder door
 - opening, 8-1, 8-2
- recorder modes, 3-10
- recorder operating controls, 1-3
- Recording, 3-9
- recording
 - automatic, 3-9
 - recording on Alarms, 3-10
 - recording on Charge, 3-9
 - recording on Mark, 3-9
 - recording on Shock, 3-9
- removing power from the defibrillator, 1-8
- replacing the battery, 8-5
- resetting the selected energy level, 4-4
- RESTORE FACTORY SETTINGS, A-8, A-9
- Review, 2-4
- review, 1-4
- reviewing Event Summary, 3-8
- R-wave, 4-5
 - quality, 4-3
- R-wave detection in cardioversion, 4-1
- R-wave marker, 4-3

- S**
- safety, 1-8
- safety symbols, iv
- SaO₂, 6-1
- selecting ECG source, 3-6
- selecting energy, 1-2, 2-1, 2-2
- SETUP LOST, 7-2, A-9
- SETUP MENU 1, A-8, A-9
- SETUP MENU 2, A-8, A-9
- setup menu 2, 3-9
- setup menus
 - configuring with, A-7
- setup mode, A-9
- setup/diagnostic menu, A-9
- shock, 2-3, 4-5
 - avoiding while pacing, 5-1
 - button, 7-10
 - event, 3-8
 - recording, 3-9
 - synchronized, 4-5
- Shock buttons, 1-2, 1-3
- shock hazard, iv
- shocking the patient, 2-3
 - recordings, 3-10
- shutdown warning, 1-8

- soap and water
 - for cleaning, 8-6
- specifications
 - synchronized shock, 4-1
- SpO₂, 1-6
 - compared to SaO₂, 6-1
 - controls, 1-6
- SpO₂ Alarm, 1-6
- SpO₂ Alarms
 - violation, 3-8
- SPO2 CABLE OFF, 7-6
- SPO2 FAILURE, 7-6
- SPO2 LIGHT INTERF, 7-6
- SPO2 LOW SIGNAL, 7-6
- SpO₂ monitoring, 6-1
- SPO2 NOISY SIGNAL, 7-6
- SpO₂ On/Off, 1-6
- SPO2 SENSOR FAIL, 7-6
- standard paddles, 2-4
- Start/Stop, 1-5
- STOP PACER, 7-5
- stored energy discharge, 4-4
- supplies, 7-9, 8-8
- switched paddles, internal, 2-6, 2-8
- switchless paddles, internal, 2-6
- SYNC, 4-2
- Sync button, 1-2, 1-3, A-8
- Sync Light, 1-7
- synchronized cardioversion, 4-1, 4-5
 - with pads, 4-2
- synchronized shock, 4-5
 - mode, 4-5
 - specifications, 4-1
- synchronized shock mode, 1-2, 1-3, 1-7
- synchronous mode, 3-5
- SYSTEM FAILURE, 7-1

- T**
- TEST 100J FAILED, 7-11
- test discharge
 - recordings, 3-10
- Time, A-10
- troubleshooting, 7-1

- U**
- Upper Alarm Limits, A-9
- USE LEADS, 4-2, 7-2
- using the defibrillator after, 4-5

- V**
- ventricular fibrillation
 - avoiding, 4-1

Reader Comment Card

We welcome your evaluation of this manual. Your comments and suggestions help us improve our publications.

Manual Title and/or Part Number _____

Please circle one number for each of the statements below.

	Strongly Agree	Agree	Neutral	Disagree	Strongly Disagree
■ The manual is well organized.	1	2	3	4	5
■ I can find the information I want.	1	2	3	4	5
■ The information in the manual is accurate.	1	2	3	4	5
■ I can easily understand the instructions and procedures.	1	2	3	4	5
■ The manual is clearly written.	1	2	3	4	5
■ The manual contains enough examples.	1	2	3	4	5
■ The examples are appropriate and helpful.	1	2	3	4	5
■ The layout and format are attractive and useful.	1	2	3	4	5
■ The illustrations are clear and helpful.	1	2	3	4	5

Please feel free to write additional comments, particularly if you disagree with a statement above. Use additional pages if you wish—the more detailed your comments, the more useful they are to us.

Comments: _____

Name: _____

Title: _____

Company: _____

Address: _____

City/State/Zip: _____

Phone: _____

Please tear out and mail in.



Fold Here
Tape edges for mailing



BUSINESS REPLY MAIL
FIRST CLASS MAIL PERMIT NO. 85 MCMINNVILLE, OREGON

POSTAGE WILL BE PAID BY ADDRESSEE

HEWLETT-PACKARD COMPANY
Cardiology Business Unit
ATTN: Quality Assurance
1700 South Baker Street
McMinnville, OR 97128-9196

NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES



M1722-91908

Copyright © 1992, 1993, 1994
Hewlett-Packard Company
Printed in USA November 1994